
www.HEBUmedical.de

HEBU HFIT®-T

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GB

HEBU HFiT®-T
Radio Frequency Thermo Ablation

> A secure and low complication rate method for inoperable tumor diseases

Relying on disclosures by the Robert-Koch-Institute, Berlin, annually more than 5.000 people in Germany are developing liver cancer. Only a few patients can be healed by a surgical intervention such as liver transplantation or liver resection. Other treatments such as chemotherapy or radiotherapy are not only complex but also are accompanied by side effects. Here out further therapeutic methods for inoperable liver cancer were developed to damage the tumor locally and conserve the surrounding tissue withal. One of the most


successful local ablative therapy methods is the Radio Frequency Thermal Ablation, so called RFA. Out of the dedicated research and development in this field a lot of valuable contributions were gained to improve the treatment of liver cancer patients. One considerable example represents the High Frequency Thermo Ablation also called HFiT®. This kind of treatment serves a minimally-invasive, secure, gentle and low-risk method for the patient to treat inoperable cancer, especially liver cancer.

Heat plus saline solution - a patented formula for success

The success of this treatment depends not only on the duration of the heat exposure but also on the electrical conductivity of the targeted tissue. To improve

conductivity, it needs to be ensured that saline solution is used for perfusion. Thus a complete ablation of the tumor cell ensues by means of the HF-Energy.

HEBU HFiT®-T is simple and superior - High Frequency Induced Thermo Therapy (HFiT®-T)

Developed and clinically tested in cooperation with German University Clinics, the HFiT® is a saline soak based further development of the RFA method. Gastroenterologists, radiologists, liver surgeons and oncologists all over the world trust this demonstrably efficient treatment option. Due to the use of saline solution heat accumulation and exsiccation of the tissue are avoided furthermore the electrical conduction and the temperature around the tip of the needle electrode are stabilized.

> Simple Handling - From needle ablation trough to HEBU HFiT®-T system application.


HEBU HFiT®-T is as flexible as the demands placed. The user is provided with a variety of needles with diverse shaft and active lengths.

The selection of the correct needle depends on the medical diagnosis, brought out of imaging methods.

HEBU HFiT®-T needles provide the following specific advantages:

- Fenestrated needle with a sharp trocar tip, similar to a biopsy needle
- Insulated needle shaft with ring markings
- Rigid electrode tip without additional moving parts

All these properties lead to a simple navigation of the needle. Risks of displacement are minimized and difficulties during the pull out of the needle out of the targeted tissue immediately after thermal ablation are eliminated.

HEBU HFiT®-T - a superior concept

During the redesign of the device, particular the focus was placed on a simple and safe handling of the device. The clearly arranged color-touch-screen enables the user to dedicate oneself with all attention to the patient.


➤ HEBU HFiT®-T for monopolar and bipolar Radio Frequency Thermo Ablation

The application technique of Radio Frequency Ablation (RFA) that can be used depends on various factors:

- Distance of the tumor to tender tissue structure
- Size and geometry of the tumor-necrosis zone to be treated.
- Active and passive components of the HF-cycle


The HEBU HFiT®-T enables the operator to use both, the monopolar and the bipolar technique.

Bipolar Thermo Ablation offers significant advantages for safe surgery in areas near heat sensitive tissues or at insufficient imaging. The neutral electrode becomes unnecessary. This system is also suited for surgeons who prefer to combine the monopolar and bipolar technique in one single session. The short distance, which does not exceed 30mm, between the tips of the needles, transfers the entire bundled HF-Energy on the tissue.



Monopolar RFA-Application

A large surface neutral electrode with low heat capacity is attached to the patient as an electrical pole. After the local or general anesthesia the thermal active electrode is inserted.

Bipolar RFA-Application

Is applied on the targeted tissue through two parallel inserted active needle electrodes. This application supersedes the need of the neutral electrode.

➤ HEBU HFiT®-T – the patient is the center of attention

Imaging became a benchmark for progress of modern medicine. HEBU HFiT®-T is compatible to most of the imaging systems (Sonography, CT). Percutaneous or intraoperative placement of the needle as well as its thermal effect can be monitored over a standard imaging equipment.

Variation of image quality may occur, depending on the imaging system. The choice of the imaging system is related on respective Indications:

- Type of therapy
- Type of anesthesia
- Experience of the operator in surgeries under image control

➤ Advance through cooperation

The HEBU HFiT®-T was enhanced in cooperation with the ultrasound department (Director: Prof. Dr. D.Strobel) of Medical Clinic 1 of the University Hospital Erlangen. The assortment of diverse Needles with

different active electrode lengths in conjunction with variable adjustable performance parameters ensure a treatment result that is optimized individually for the patient. The ultrasound images are showing a

direct comparison before and after a sonographic controlled RFA on a patient with hepatocellular carcinoma (HCC) in the right hepatic lobe.

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HF-Generator

HEBU HFiT®-T

Art.-No.	HF 9518
Power Supply	100 – 260 Vac, 50-60 Hz
Power Consumption max.	108 VA (100 Vac) / 100 VA (230 Vac)
HF-Power output	60 W (150 Ohm) 268 Vss/95 Veff
Nominal HF-frequency	500 kHz
Adjustable Parameters	HF-Power in W, Treatment time in s, Flow rate in ml/h, Temperature in °C (only when using the temperature probe), Synchronization (P-Z-F) ON/OFF, Automatic (P-F) ON/OFF, Program memory for nine Programs.
Visualization of applied Parameters	Real-Time Impedance monitoring, Total Energy input, Operating temperature (only when using the temperature probe), Voltage in V, Current in mA, Power in W, Impedance in Ohm
Dimensions	L x H x W = 350 x 370 x 145 mm
Weight	8,1 kg

Perfusion pump

Art.-No.	HF 9518-01
Power Consumption max.	23 VA
Feed rate	0,1 – 200 ml/h for 50ml syringe, controlled over HF-Generator via RS 232 interface
Dimensions	L x H x W = 330 x 120 x 155 mm
Weight	2,4 kg

Trolley

Art.-No.	HF 9529
Dimensions	L x H x W = 420 x 420 x 850 mm
Weight	28 kg

Norm

EU Directive 93/42/EEC (CE 0123)

Ordering Data

HEBU HFiT®-T-System

HF 9518-00	HEBU HFiT®-T with Temperature control (HF 9518), Perfusion Pump (HF 9518-01), RS 232 Interface cable (HF 9418-02), Trolley (HF 9529), Double foot switch for HEBU HFiT®-T (HF 9537-02)
------------	--

Needle Electrodes, sterile packed

HF 9945-15	HEBU HFiT®-T Needle Electrode, \varnothing x shaft length x electrode length: 2,0 x 150 x 15 mm
HF 9945-20	HEBU HFiT®-T Needle Electrode, \varnothing x shaft length x electrode length: 2,0 x 150 x 20 mm
HF 9946-15	HEBU HFiT®-T Needle Electrode, \varnothing x shaft length x electrode length: 2,0 x 200 x 15 mm
HF 9946-20	HEBU HFiT®-T Needle Electrode, \varnothing x shaft length x electrode length: 2,0 x 200 x 20 mm
HF 9946-25	HEBU HFiT®-T Needle Electrode, \varnothing x shaft length x electrode length: 2,0 x 200 x 25 mm
HF 9946-30	HEBU HFiT®-T Needle Electrode, \varnothing x shaft length x electrode length: 2,0 x 200 x 30 mm

Accessories

HF 9518-01	Perfusion Pump for 50 ml syringe, 0,1 – 200 ml/h
HF 9418-02	RS 232 Interface cable
HF 9529	Trolley
HF 9537-02	Double foot switch for HEBU HFiT®-T-system, class AP, 6 Pin plug, cable length 4,5 m
HF 9564 H	Single-use security electrode neutral with double separated contact surface, Hydrogel, Conductive area 110 cm ² , 50 pcs per carton
HF 9949-00	Temperature probe, Type K, \varnothing 0,5 mm, with cable, 3 m, sterilisable
HF 9949-01	Y-connector for insertion of the Temperature probe into the HEBU HFiT®-T Needle and for connection of the perfusion hose, 5 pcs per carton
HF 9949-02	Connection hose with Luer-Lock-Connection, 1,9 m, 25 pcs per carton
HF 9949-03	Syringe 50 ml, 50 pcs per carton

Cables

HF 9526-05	Potential equalization cable, 0,5 m
HF 9526-10	Potential equalization cable, 1,0 m
HF 9526-30	Potential equalization cable, 3,0 m
HF 9569-30	Cable for disposable return plates, 3,0 m
HF 9569-50	Cable for disposable return plates, 5,0 m
HF 9571-30	Monopolar cable for HEBU HFiT®-T Needles, sterilisable, 3,0 m
HF 9572-30	Bipolar cable for HEBU HFiT®-T Needles, sterilisable, 3,0 m
HF 9979-01	Mains cable, e.g. for Germany, France, Austria, Scandinavia, 2,5 m
HF 9979-02	Mains cable for Switzerland, 2,5 m
HF 9979-03	Mains cable, e.g. for Taiwan, Korea, USA, Canada, 2,5 m
HF 9979-04	Mains cable for Italy, 2,5 m
HF 9979-05	Mains cable for Great Britain, 2,5 m
HF 9979-06	Mains cable for Argentina and Australia, 2,5 m

HEBUmedical GmbH
Badstraße 8
78532 Tuttlingen | Germany

Tel. +49 7461 9471 - 0
Fax +49 7461 9471 - 22

info@HEBUmedical.de
www.HEBUmedical.de