HF- Zubehör / HF accessories / accessoires HF / accessorios AF / Accessori RF

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Seite 3

Operating Manual

Page 13

Mode d'emploi

Page 21

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Página 31

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Página 41

DEUTSCH

ENGLISH

FRANÇAIS

ESPAÑOL

ITALIANO

Inhalt

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Contents

1	Used Symbols	13
2	Scope	13
2.1	Inteded Use	13
2.2	Contraindication	15
3	Storage	15
4	Cleaning and sterilization	15
5	Instructions for Use	16
5.1	Connection to HF unit in monopolar mode.....	16
5.2	Connection to HF unit in bipolar mode	17
6	Warnings	17
7	Safety Remarks	18
8	Additional information for the use of argon accessories	18
8.1	Intended use.....	18
8.2	Instrument preparation	18
8.3	Connecting the argon surgical handle.....	19
8.4	Safety information	19
8.5	Open surgical application	19
8.6	Laparoscopic application.....	20
9	Servicing and manufacturer address	20

1 Used Symbols

Symbol	Definition
	CE-marking
	Attention
	Keep away from sunlight
	Store at dry place
	No intervention on patients with cardiac pacemakers.
	Manufacturer

2 Scope

2.1 Inteded Use

The HF-accessories include the following class I and IIb products:

- **HF-cables**
A HF-cable is a product that acts as electrical conductor between an electrosurgical unit and an electrode, which is used for transmission of energy. The energy may contain information's or not.
Maximum electrical capacity: 7,5kV_p
- **Handles**
Handles are used to connect the assistance instruments for monopolar or bipolar coagulation on HF-units.
Maximum electrical capacity: 7,5kV_p
- **Electrode handles without pushbuttons:** Handles are used to connect auxiliary instruments for monopolar or bipolar coagulation to HF devices. It is a reusable article.
Electrode extensions; Adaptor for hand switch pencils: Adapters and HF extender are usually small pieces that are used to connect the lead of the electrodes to an HF device. They are reusable items.
HF Cables: An HF cable is a product that acts as an electrical conductor between an electrosurgical unit and electrode and is used for the purpose of transmitting energy. This transfer of energy may or may not contain information. It is a reusable product.
- **HF-Handles with pushbuttons:** Pushbuttons are used to connect auxiliary instruments for monopolar or bipolar coagulation to HF devices and are intended as a contact and switching element for HF electrodes. Activation takes place by means of a button on the HF handle.

They are reusable items.

Foot switches for electrosurgical unit: A footswitch (pedal) is a switching element for monopolar and bipolar coagulation on HF devices. The activation takes place by means of pedal. It is usually operated by the surgeon or doctor to address the HF device to which it is connected. It is a reusable article.

- Elektroden

Monopolar electrodes:

- Disposable monopolar electrodes, Ø2,4mm
- Monopolar electrodes for dental and cosmetic applications, Ø1,6mm
- Monopolar standard electrodes, Ø2,4mm/ Ø4mm
- Very delicate, flexible electrodes for the plastic surgery, Ø2,4mm
- Monopolar electrodes with insulated shaft, Ø2,4mm/ Ø4mm
- Microdissection needle electrodes, 2,4mm/ Ø4mm
- Needle electrodes with ceramic isolation, 2,4mm/ Ø4mm
- Monopolar electrodes for arthroscopy, Ø2,4mm/ Ø4mm
- Electrodes for laparoscopic surgery, with heat resistant ceramic tip isolation
- Monopolar Argon electrodes, Ø4mm

Maximum electrical capacity: 4,3kV_p

Bipolar electrodes:

- Bipolar coagulation electrodes for ENT
- Bipolar coagulation electrodes for Laparoscopy, Ø3mm/ Ø5mm
- Bipolar arthroscopy electrode

Maximum electrical capacity: 1,0kV_p

- Monopolar und bipolar Instruments

Monopolar and bipolar Instruments are devices covered with a sheath of electrical insulation, except at the distal working tips, designed for use in open surgery for the coagulation of vessels or tissue bundles to prevent bleeding. They are connected to either an electrosurgical diathermy or an electrocautery generator which provides the electrical current to the device. They are typically constructed as: 1) long, tweezers-like instruments with tapered, flat, blades which are closed by the surgeon applying pressure to the conjoined legs; or 2) scissors-like instruments with long, curved, clamp-like blades. They are made of metal and synthetic materials in a variety of shapes and sizes. These are reusable devices.

Monopolar Instruments:

- Instruments for polypectomy
- Instruments for papillotomy
- Monopolar Instrument with exchangeable electrodes
- Monopolar forceps with and without socket
- Flexible forceps for hot biopsy

Maximum electrical capacity: 4,3kV_p

Bipolar Instruments:

- Bipolar coagulation forceps for laparoscopy
- BiSeal® Vessel sealing instrument
- Bipolar forceps, insulated with exchangeable tips
- Bipolar suction forceps, insulated
- Bipolar scissors, insulated
- Bipolar clamps

Maximum electrical capacity: 0,5kV_p

- **HF coagulation suction tube**

A Monopolar suction tube for suction of fluids and for coagulating soft tissue during surgery.
Maximum electrical capacity: 4,3kV_p

- **Adaptors for HF-handles, HF- extensions**

A connecting device, usually a small piece used to join the cable from electrodes (active or return) to an electrosurgical unit.
Maximum electrical capacity: 7,5kV_p

When used as intended, these products are more or less subject to wear, depending on the intensity of use. This kind of wear is caused by technical factors and cannot be avoided. Therefore carry out a thorough visual examination and a function test prior to each application of these products. If a new product has externally visible defects or does not operate according to the instructions, contact the manufacturer or his responsible representative immediately.

2.2 Contraindication

In the following cases, no or limited HF surgery should be performed:

- Patients with implanted pacemakers or pacemaker leads

In areas where flammable or explosive agents are present, e.g. Skin cleansing, degreasing and disinfecting agents but also endogenous gases, e.g. in the gastrointestinal tract (risk of fire and explosion)

- Acute emergency situation
- serious coagulation disorders
- Severe impairment of lungs or cardiovascular system

The HF neutral electrodes must not be used on damaged and injured skin.

Die HF-Monopolare Pinzetten, HF-Bipolare Pinzetten, HF-Monopolare Klemmen, HF-Bipolare Klemmen, HF-Bipolare Elektroden und HF-Monopolare Elektroden

The HF bipolar clamps, HF Monopolar electrodes, HF bipolar electrodes, HF monopolar forceps, HF bipolar forceps and HF monopolar clamps should not be used in conjunction with the central nervous or circulatory system. Tube ligation is not a safe application.

HF bipolar scissors should not be used in patients with heart-cutters.

3 Storage

Temperature: -20 °C - +50°C

Relative humidity: 0 -75%, not condensing

Atmospheric pressure: 500 -1600 mbar

Products must be protected from direct sunlight in storage. The storage of the products in their original packaging until their first use is recommended.

4 Cleaning and sterilization

Products are supplied in an unsterile state and must be cleaned and sterilised before use.

The disposable monopolar electrodes and electrode handles are sterile packed and must not be cleaned and sterilized before use. They must be disposed of after single use.

Do not clean HF cables, adapters or electrode handles in an ultrasonic bath.

Do not sterilise HF cables, adapters or electrode handles and electrodes in hot air sterilisers.

If a Washing-Desinfection unit is used for cleaning: max 95 °C. The HF cable or handle also may be immersed in a combined cleaning and disinfection solution (e.g. neodisher® MediClean forte, Dr. Weigert Co.). In such a case, the service life may then be reduced due to oxidation of the contacts and crystallisation in the plugs. Please observe the chemical manufacturer's instructions. Then rinse cable, connectors, and handle or electrode thoroughly with demineralized water. Steam-sterilize products at 121°C / 20 min minimum dwell time or 132 resp. 134°C / 4min minimum dwell time (max. 10 min).

5 Instructions for Use

Thoroughly examine all products for visible irregularities. Do not use cables with a brittle or defective insulation. In case of doubt check operation by means of a suitable testing device.

5.1 Connection to HF unit in monopolar mode

First fix the neutral electrode ("return plate") to the patient, preferably on the upper arm or on the thigh. Ensure that the skin is free from hair and grease. If a reusable neutral electrode (return plate) is used, do not apply conductive gel. The current paths inside the patient's body should be short and must run diagonally. Current paths must never run in a transverse direction through the body or across the thorax. The patient must be insulated against all electrically conductive parts. Ground the operating table. Place the patient on a dry electrically insulating layer. Prevent skin areas touching each other. Place dry gauze in between. Then connect the neutral electrode, the foot switch and the "active" HF cable to the HF unit. In order to avoid interference with camera system, do not run HF cables directly in line with camera cables. Never place HF cables directly on the patient's skin. Never lay HF cables in loops. Next, connect the instrument or handle with coagulation electrode, and switch on the HF unit. The electrode with the proper shaft diameter (4 mm or 2.4 mm) must be firmly attached to the handle tip. For handles with finger switches:

Yellow switch = actuation of cutting current

Blue switch = actuation of tissue coagulation

Always refer to the instruction manual of the HF unit and to the general information it provides relating to electrosurgical procedures.

5.2 Connection to HF unit in bipolar mode

In bipolar HF surgery the coagulation is restricted to the area between the two electrodes of the bipolar instrument. Consequently, a neutral electrode is necessary. Otherwise, follow the recommendations regarding patient preparation in monopolar operating mode and the corresponding instructions for the HF unit. When using the bipolar “Auto start” mode, the cable used should not be longer than 4.5 m. Alternatively, use a twin cable with max. length of 5m.

6 Warnings

The products may only be used in combination with HEBU accessories and only by clinically instructed and qualified personnel. Electrodes must not be twisted, as this may cause damage to the material and the electrode might break. Power output of the HF unit: Always select the lowest possible output of the HF unit for any procedure. If the coagulation capacity of the electrode is below normal, do not increase the output of the HF unit without a thorough preliminary check. In particular, the following must be checked: correct attachment of all HF cables and plugs, correct activation of current by the finger - or footswitches, undamaged insulation of the HF cables and instrument, and check that the distal end of the electrode is dean and undamaged. Before operating, check that the signal transmission of the monitors works without interference (e.g. without noise); this check is made by pressing the yellow button (CUT) and the blue button (COAG).

Always unplug the HF-cable by holding on the connector. **Never pull the cable itself**, as this may lead to defects. Cables with defective/brittle insulation or with a broken electrical wire (due to severe bending) might cause a fire. Active HF instruments and HF handles must not be placed on the patient in order to avoid injuries in case of involuntary HF activation or activated “Auto Start” mode.

For cleaning the active electrode, you must deactivate the “Auto Start” mode (only in bipolar applications) or disconnect the instrument from the generator. Directly after activation the tips could be hot. This may lead to burns.

Endogenous burns are those caused by high current density in the patient’s tissue. Possible reasons are: The patient is inadvertently positioned in contact with electrically conductive parts. Direct contact between skin areas and HF cables may lead to electrical capacitance which in turn may cause burns.

Exogenous burns are those caused by heat of ignited fluids or gases. They may also be caused by explosions. Possible reasons include: Ignition of skin cleaning agents and disinfectants, ignition or narcotic gases etc.

Cardiac pacemakers may be damaged by electrosurgical current. Consult a cardiologist before operating. Never use electrosurgical current on outpatients fitted with pacemakers.

7 Safety Remarks

!	Check the isolation on the electrode and the shaft extensions for damage before each use.
!	Do not rebend needles, loops or knives. Spatulas can be bent.
!	If damaged, do not use these products!
!	Do not use in the presence of combustible or explosive materials!
!	Never lay instruments on the patient or in his/her direct vicinity!
!	HEBUmedical expressly warns against modifying the product. Any modification exempts HEBUmedical in from any and all liability. .
!	Protect these products from any form of mechanical damage! Do not throw! Do not: use force!
!	Defective HF cables or electrode handles may not be repaired. Always replace a defective part with a new one.

8 Additional information for the use of argon accessories

8.1 Intended use

The argon surgical handle is a re-usable electrosurgical handle with gas feed. It is used in argon gas assisted HF surgery. It is suitable for the fitting of monopolar, rigid argon electrodes, which are used for argon assisted coagulation or for argon gas assisted electrosurgical cutting. The cutting and coagulation currents can be controlled on the HF surgical unit by means of two finger switches. Fundamentally, the argon surgical handle can also be operated by means of a foot switch. These products may only be used in medical establishments by trained personnel. Read the instruction manual carefully before using the products.

8.2 Instrument preparation

Mechanical precleaning of argon knife and needle electrodes can be carried out in the same way as with conventional surgical electrodes. As liquid can enter the gas feed pipe of the handle during the cleaning process, the lumen of the gas feed pipe should be blown free with compressed air, max. pressure 4 bar, after cleaning. For this purpose, the compressed air fitting is placed loosely on the Luer Lock connector of the argon surgical handle (do not blow out from distal to proximal). If a Washing-Desinfection unit is used for cleaning: max 95 °C. The argon surgical handle is sterilized in the autoclave at 121°C / 20 min minimum dwell time or accordingly 132 resp. 134°C / 4min minimum dwell time (max. 10 min).

Do not sterilize in hot air!

Gas sterilization is not recommended, as the airing times have not been evaluated for this product.

8.3 Connecting the argon surgical handle

It must be ensured first of all that the HF generator and the argon gas control unit are ready for operation. The gas bottle valve must be open and there must be an adequate supply of gas.

1. Remove the handle and electrodes from the container in a sterile manner.
2. Visually inspect the instrument for damage to the insulation. Do not use damaged instruments! HEBU warns against repairing damaged areas. The instrument is operated at high coagulation voltages!
3. Connect 3-pin plug to the HF surgical unit.
4. Plug the required argon electrode into the electrode socket on the handle. Make sure, that the plastic hexagonal section of the electrode shaft is fully inserted into the handle.
5. Make the argon gas connection to the argon control unit. To do this, connect the male Luer-Lock connector of the handle to the female Luer-Lock gas outlet of the argon control unit.
6. Activate the purge function of the argon gas control unit in order to flood the gas feed pipe of the handle with argon gas.
7. Operate the foot switch and/or finger switches of the argon surgical handle one after the other and check the operation of the system (activation of the argon gas control unit and HF generator).

8.4 Safety information

Rigid argon coagulation electrodes are operated exclusively with spray coagulation current.

Always operate argon coagulation electrodes by means of the blue buttons!

Basically, argon knife and needle electrodes (cutting electrodes) can be used with argon gas assistance with all type of monopolar cutting and spray coagulation current. Unmodulated cutting current is recommended for gas-assisted cutting.

In order to “fire” the argon plasma, the argon electrode is moved within ca. 1 cm of the tissue. During the argon-assisted coagulation process (beaming), this distance can be increased to approx. 2 – 2.5 cm before the plasma extinguishes.

8.5 Open surgical application

The argon surgical handle allows a gas feed in the range 0-9.5 l/min. The gas flow is selected on argon control unit depending on the selected argon electrode and the indication. Flow rates of 2-4 l/min have been shown to be beneficial for argon electrodes in open surgical applications. In this case, the output power of the HF generator in spray mode is ca. 40 W.

When cutting an output power of approx. 60 W has been shown to be beneficial for argon needle electrodes and a power output of approx. 90 W for argon knife electrodes. In the cutting mode, the gas flow is ca. 2-4 l/min.

8.6 Laparoscopic application

The trocar valve must be open during argon application in order to ensure that a critical intra-abdominal pressure is not exceeded. If an insufflator pressure has fallen below the critical value once more.

Whenever the intra-abdominal pressure is monitored by an insufflator, an independent, continuous monitoring of the intra-abdominal pressure by the operating team is required. The use of a pressure-relief controlled insufflator is recommended.

9 Servicing and manufacturer address

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HEBUmedical GmbH