HF 9928-01 HF 9928-02
HF 9928-03 HF 9928-04
HF 9928-05 HF 9928-06
HF 9929-01 HF 9928-07
HF 9929-02
HF 9929-03
HF 9929-04

APC-Sonde / APC probe / Sonde APC / Sonda APC / Sonda APC

**Wiederverwendbar / Reusable / Réutilisable/ Reutilizable /
Riutilizzabili**

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Seite 3

Operating Manual

Page 13

Mode d'emploi

Page 23

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Página 33

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Página 43

DEUTSCH

ENGLISH

FRANÇAIS

ESPAÑOL

ITALIANO

HF 9928-01 *HF 9928-02*
HF 9928-03 *HF 9928-04*
HF 9928-05 *HF 9928-06*
HF 9929-01
HF 9929-02
HF 9929-03

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Flexible Probe for Argon Plasma Coagulation- Reusable

Contents

1	Symbols	14
2	Scope	14
2.1	Inteded Use	15
2.2	Contraindications	15
3	Storage.....	15
4	Preparation.....	16
4.1	Cleaning	16
4.2	Sterilization	17
5	Connecting the probe	18
6	Safety remarks	18
7	Application remarks	20
8	Manufacturer and servicing address	20

1 Symbols

Symbol	Definition
	CE-marking
	Lot number
	Order number
	Attention
	Note instructions for use
	Temperature limitation
	Store at dry place
	Keep away from sunlight
	Manufacturer
	Non-steril

2 Scope

The use of reusable APC problems in the following areas:

- General medicine
- Gynecology
- liver surgery
- Traumatology
- Cardiac and Thoracic Surgery
- ENT

Example of products:

HEBU APC probe [1] HF 9928-01

Connecting cable [2] HF 9928-00

Rinsing adapter [3] HF 9929-98

2.1 Inteded Use

HF Argon coagulation probe; HF-HEBU APC-Probe: The flexible probes for argon plasma coagulation (APC probes) are used for argon gas-assisted coagulation and are used in oesophago-, gastro-, broncho-, colono- and rectoscopically assisted interventions in the gastrointestinal tract or in the tracheo-bronchial space. The APC probes can be used in conjunction with electrosurgical units (high frequency generators) and an argon unit with an output voltage of max. 4 kVs are operated. They are connected to the unit with the reusable and steam sterilizable cable. The probes are reusable and steam sterilizable. The activation of the gas supply as well as the HF-current takes place via the footswitch of the high frequency generator.

HF Argon Probe – rinsing adapter: A connector for cleaning an HF argon coagulation probe. The rinsing adapter has a Luer Lock connection. It is a reusable product.

2.2 Contraindications

There are no known contraindications.

3 Storage

Store the flexible APC probes and their connecting cable in a cool and dry location. Protect from the effects of sunlight.

Medical products should never be stored or operated in dusty environments on principle.

Storage temperatures: -25° to +70°C

Relative humidity: 10 to 75 %
Air pressure: 500 to 1060 hPa

4 Preparation

**Instruments are not sterile in their as-delivered condition!
Clean and sterilize the product prior to first application!**

Remark:

Observe national statutory regulations, national and international standards and directives and the user's own hygiene regulations for the preparation of instruments.

Note that successful preparation of this medical product can only be ensured following prior validation of the preparation process. This is the responsibility of the operator / preparatory staff member.

4.1 Cleaning

Before cleaning the APC probe, disconnect from the connecting cable.

Prior treatment

The APC probe must be rinsed through with water after using on a patient (drinking water quality or sterilized water). Should this not be available in exceptional circumstance, the probe must at least be blown through using an air-filled syringe. Ensure that no kinks are formed in the probe's tube.

The probe must not be rinsed through with disinfectant solutions, as residues of this can impair the cleaning result as a result of interaction with the treatment medium during the mechanical preparation process.

Laying the instrument in cleaning or disinfectant solution can severely reduce its service life.

Note! If the described prior treatment is not carried out, the success of subsequent mechanical cleaning processes cannot be guaranteed.

Mechanical cleaning

HEBU APC probes and their connecting cable can be mechanically prepared. Thermal disinfection at 93°C is possible. For this purpose, connect the HEBU APC probes to nozzles using a special machine insert (such as an injector trolley, MIC insert) in order to ensure that they are **rinsed through during the cleaning process**. Programs which perform thermal disinfection during the final rinse cycle are suitable for mechanical preparation.

If there is no suitable hose adapter available in the rinsing machine, the supplied rinsing adapter with Luer lock connector is available.

A suggested program encompassing the necessary main steps is outlined below:

Program step	Water requirement quality	Treatment medium
Pre-rinse	At least drinking water	None
Cleaning	At least drinking water	Mildly alcoholic cleaning agent
Neutralization	At least drinking water	Neutralization component
Intermediate rinse 1	At least drinking water	None
Intermediate rinse 2	At least drinking water	None
Thermal disinfection / final rinse	Fully demineralized water	None

Treatment medium

For cleaning, pH neutral and mildly alkaline cleaning agents approved for the preparation of surgical instruments can be used. For neutralization, citric acid or phosphoric acid-based agents can be used.

Cleaning temperatures

The temperature should **not exceed 50-55°C** during the cleaning phase.

The information from the cleaning agent manufacturer and the manufacturer of the automatic cleaning machine must be observed.

Avoid contact with sharp objects!

As fluid can penetrate into the probe's tube during cleaning, after the cleaning program has terminated blast through the lumen using compressed air. To do this, loosely place the compressed air attachment on the probe's plug (do not blast through from distal to proximal).

4.2 Sterilization

The HEBU APC probe and the connecting cable are sterilized at a temperature of up to +134°C.

Do not sterilize the probe and connecting cable in an assembled condition!

The probe may not be sterilized in hot air!

Gas sterilization is not recommended, as the airing times have not been evaluated for this product.

5 Connecting the probe

Initially, ensure that the electrical surgical device and the (integrated) argon unit are ready for operation. Ensure that the gas cylinder valve is open and that there is an adequate supply of gas (> 30 bar).

1. Remove the connecting cable (p. 3, [2]) and HEBU APC probe (page 3, [1]) from the sterilization container under sterile conditions.
2. Carry out a visual inspection of the instrument for possible damage to the insulation. Do not use damaged instruments! HEBUmedical warns users not to attempt to repair or patch up damaged areas. The instrument is operated at high coagulation voltage levels!
3. Check whether the working channel of the endoscope is suitable for the required probe diameter (3.2 mm, 2.3 mm or 1.5 mm).
4. Connect the three-pin plug of the connecting cable to the HEBU electrical surgery device.
5. Insert the required probe with plug into the coupling at the end of the connecting cable. Ensure that the plug is pushed into the coupling as far as it will go.
6. Establish the argon gas connection to the electrical surgery device / the argon unit. To do this, connect the female Luer lock connector of the connecting cable to the male Luer lock gas outlet of the electrical surgery device / the argon beamer. If required, an intermediate sterile filter can be connected.
7. Activate the function “C” at the electrical surgery device / the argon unit in order to flood the gas lumina of the connecting cable and the probe with argon gas.
8. Actuate the blue foot pedal and check the system function (activate the argon beamer and generator).

The probe is ready for operation.

6 Safety remarks

! Never place the gas outlet nozzle of the HEBU APC probe directly on the tissue during gas activation! **Danger of gas emphysema!**

! Do not guide the gas outlet nozzle vertically over the tissue during activation. An angle of around 45° to the tissue has proven successful.

! Do not blast argon into the vessel system. **Danger of embolism!**

	<p>Always hold the applicator end within the visual range of the endoscope during activation! Never activate the coagulation stream without visual inspection!</p>
	<p>Endogenous gas in the gastrointestinal tract can be combustible or explosive. Before any electrosurgical application, clear any such gases in suspected areas by flushing with argon gas (e.g. by activating the beamer on the HF output setting "0" or using the purge function "C/Argon Flow" of the beamer).</p>
	<p>Argon plasma can ignite combustible materials (e.g. the plastics used in the endoscope or tubes). This can result in a hazard, particularly in oxygen-enriched environments!</p>
	<p>The distal end of the probe must project by at least 10 mm from the working channel of the endoscope. There are marking rings provided at the distal probe end for orientation purposes.</p>
	<p>Do not introduce oxygen or other combustible gases or fluids before and during application in the tracheo-bronchial system. During protracted beamer applications, carry out alternating ventilation!</p>
	<p>In the event of laparoscopic applications, the argon gas flow brings about an increase in intra-abdominal pressure. The use of electronic pressure release-controlled CO₂ insufflators is recommended. During activation, monitor intra-abdominal pressure and where necessary release pressure through the trocar!</p>
	<p>For endoscopic applications with video endoscopes, do not direct the argon gas jet towards the camera chip.</p>
	<p>Protect the probe from mechanical overloading. Ensure that the hose does not kink!</p>
	<p>As high HF voltage levels are required for gas ionization, pay particular attention to ensuring that the electrical insulation of cables and applicators is undamaged.</p>
	<p>Electrode wires which project out from the flexible tube at the distal end can result in mechanical damage to tissue</p>
	<p>As argon gas-supported application is a monopolar high-frequency form of surgical application, the safety rules for monopolar high-frequency surgery must be observed!</p>
	<p>Due to the manufacturing method used, quality 4.8 argon in gas cylinders is germ free. In order to guarantee freedom from germs, the cylinder valve and the gas outlet at the gas supply device must be disinfected prior to connection.</p>

7 Application remarks

- HEBU APC probes are generally operated using a spray coagulation current mode. Select this current mode and always operate using the blue pedal!
- To "ignite" the argon plasma, approach the tissue with the distal probe opening to a distance of around 1-5mm.
During the argon-supported coagulation process, the distance can be increased to appr. 2 - 2.5 mm before the plasma is extinguished.

Flexible endoscopy and laparoscopy applications:

The diameter of the endoscope working channel should be at least 0.3 mm bigger than the specified probe diameter.

The maximum flow rate of 2 l/min must not be exceeded!

During endoscopic application, flow rates of 0.3 to 0.8 l/min are recommended at a spray output rating of appr. 30 W.

HEBU article no.	Diameter	Length	Flow rate
HF 9928-01 HF 9929-01	1,5 mm	1,5 m	0,3 - 0,4 l/min
HF 9928-03	2,3 mm	2,2 m	0,4 - 0,6 l/min
HF 9928-02 HF 9929-02	2,3 mm	1,0 m	0,4 - 0,6 l/min
HF 9928-04 HF 9929-03	3,2 mm	2,2 m	0,5 - 0,8 l/min

For the

individual probe sizes, this means:

In order to ensure that a critical intra-abdominal pressure level is not exceeded, the trocar valve must be opened during laparoscopic applications. If an insufflator pressure alarm gives off a signal, the argon application must be stopped until the intra-abdominal pressure has dropped back to below the critical level.

Whichever method is used to monitor intra-abdominal pressure by an insufflator, independent and continuous monitoring of intra-abdominal pressure by the OP team is a requirement. The use of a pressure release-controlled insufflator is recommended.

8 Manufacturer and servicing address

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HEBUmedical GmbH