HF 9945-xx
HF 9946-xx

HEBU HFit[®]-T

Nadelelektroden / Needle Electrodes/ Elettrodi ad ago

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Seite 4

Operating Manual

Page 9

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Page 15

DEUTSCH

ENGLISH

ITALIANO

HF 9945-xx
HF 9946-xx

INHALT

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HF 9945-xx
HF 9946-xx

CONTENT

1	Symbols.....	10
2	Scope.....	11
2.1	Inteded Use.....	11
2.2	Contraindication	11
3	Safety Remarks	11
4	Instructions for Use	12
5	Connection of the Needle Electrode	12
6	Types of Needle Electrodes	12
7	Storage.....	13
8	After-Sales Service Address	13

1 Symbols

Symbol	Definition
	CE-marking
	Lot number
	Order number
	Attention
	Note instructions for use
	Store at dry conditions
	Keep away from sunlight
	Manufacturer
	Do not reuse
	Sterile

2 Scope

HEBU HFiT®-T Needle Electrode is used in HF surgery for coagulation and tissue cutting.

2.1 Inteded Use

HFiT®-T needle electrodes serve as an applicator for radiofrequency thermal ablation of inoperable tumors of the liver and other pathogenic biological tissues. High-frequency current is introduced into the tissue via the HFiT®-T needle electrode and perfused with physiological saline solution. HFiT®-T needle electrodes must only be used with the HEBU HFiT®-T electro surgery system. It is a disposable product.

2.2 Contraindication

The HF monopolar electrodes must not be used in conjunction with the central nervous or circulatory system. Patients with pacemakers and in vivo cardiac defibrillators may have restrictions on use. HF monopolar electrodes should not be used in the immediate vicinity of bone, tooth root and metal.

3 Safety Remarks

	HFiT®-T needle electrodes are intended for single use only!
	Do not use HFiT®-T needle electrodes when packaging is damaged!
	Do not load HFiT®-T needle electrodes by bending or twisting!
	Do not use device in explosive environments!
	Verify that combustible materials used in cleaning agents, disinfectants or solvents for adhesives have been completely removed prior to surgery.
	Please refer to the manual instruction of the HEBU HFiT®-T to obtain information about radio frequency thermo ablation.

4 Instructions for Use

- Before inserting HFiT[®]-T needle electrode into the target area, aeration of the infusion channel and flushing with physiological saline solution must be ensured, in order to prevent blood rising into the inside of the needle.
- Air bubble inclusions must be avoided.
- The insertion can take place under ultrasound, MRT or CT monitoring.
- If the need arises, any blockage caused by clotting inside the needle can be eliminated with the sterile obturator. If this is not successful, HFiT[®]-T needle electrode must be rejected.
- HFiT[®]-T needle electrodes are supplied sterile.

5 Connection of the Needle Electrode

- The HFiT[®]-T needle electrode is connected to a standard commercial syringe pump infusion line via the Luer-Lock connection.
- The HFiT[®]-T needle electrode is contacted via the lateral socket with cable HF 9571-30.

6 Types of Needle Electrodes

Different types of HFiT[®]-T needle electrodes are available:

HEBU Item No.	Diameter	Length	Active Electrode
HF 9945-15	2,0 mm	150 mm	15 mm
HF 9945-20	2,0 mm	150 mm	20 mm
HF 9946-15	2,0 mm	200 mm	15 mm
HF 9946-20	2,0 mm	200 mm	20 mm
HF 9946-25	2,0 mm	200 mm	25 mm
HF 9946-30	2,0 mm	200 mm	30 mm

7 Storage

HFiT®-T Needle Electrodes must be stored in dry, dust free conditions and kept away from direct sun rays. Avoid temperature variation.

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HEBUMedical GmbH
HEBUmedical GmbH