# HF-Kabel, HF-Handgriffe ohne Taster

HF cables, HF handles without switch

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Seite 3

**Operating Manual**

Page 9

DEUTSCH

ENGLISH

# Inhalt

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# Contents

<b>1</b>	<b>Used Symbols</b> .....	<b>10</b>
<b>2</b>	<b>Scope</b> .....	<b>10</b>
2.1	Inteded Use .....	10
2.2	Contraindication .....	11
<b>3</b>	<b>Storage</b> .....	<b>11</b>
<b>4</b>	<b>Cleaning and sterilization</b> .....	<b>12</b>
<b>5</b>	<b>Instructions for Use</b> .....	<b>12</b>
5.1	Connection to HF unit in monopolar mode.....	12
5.2	Connection to HF unit in bipolar mode .....	12
<b>6</b>	<b>Warnings</b> .....	<b>13</b>
<b>7</b>	<b>Safety Remarks</b> .....	<b>13</b>
<b>8</b>	<b>Servicing and manufacturer address</b> .....	<b>14</b>

## 1 Used Symbols

Symbol	Definition
	CE-marking
	Attention
	Keep away from sunlight
	Store at dry place
	No intervention on patients with cardiac pacemakers.
	Manufacturer

## 2 Scope

Typical procedures with HF surgery are:

- Dermatology: epilation, wart removal
- ENT: polyp removal, monopolar and bipolar hemostasis
- Gynecology: electroconversion of the portio, sterilization by bipolar coagulation of the fallopian tube
- Urology: Transurethral resection of the prostate, bladder and bladder wall papyloomas, cutting and coagulation in transabdominal procedures
- Laparoscopy and cystoscopy: removal of polyps
- General Surgery: Primary sections, tissue sectioning, plastic surgery, hemostasis, coagulation

### 2.1 Inteded Use

The HF-accessories include the following class I and IIb products:

#### HF-cables

A HF-cable is a product that acts as electrical conductor between an electrosurgical unit and an electrode, which is used for transmission of energy. The energy may contain information's or not.

Maximum electrical capacity: 7,5kV<sub>p</sub>

#### Handles

Handles are used to connect the assistance instruments for monopolar or bipolar coagulation on HF-units.

Maximum electrical capacity: 7,5kV<sub>p</sub>

#### Adaptors for HF-handles, HF- extensions

A connecting device, usually a small piece used to join the cable from electrodes (active or return) to an electrosurgical unit.

Maximum electrical capacity: 7,5kV<sub>p</sub>

**Electrode handles without pushbuttons:**

Handles are used to connect auxiliary instruments for monopolar or bipolar coagulation to HF devices. It is a reusable article.

**Electrode extensions;** Adaptor for hand switch pencils: Adapters and HF extender are usually small pieces that are used to connect the lead of the electrodes to an HF device. They are reusable items.

**HF Cables:** An HF cable is a product that acts as an electrical conductor between an electrosurgical unit and electrode and is used for the purpose of transmitting energy. This transfer of energy may or may not contain information. It is a reusable product.

When used as intended, these products are more or less subject to wear, depending on the intensity of use. This kind of wear is caused by technical factors and cannot be avoided. Therefore carry out a thorough visual examination and a function test prior to each application of these products. If a new product has externally visible defects or does not operate according to the instructions, contact the manufacturer or his responsible representative immediately.

## 2.2 Contraindication

In the following cases, no or limited HF surgery should be performed:

- Patients with implanted pacemakers or pacemaker leads

In areas where flammable or explosive agents are present, e.g. Skin cleansing, degreasing and disinfecting agents but also endogenous gases, e.g. in the gastrointestinal tract (risk of fire and explosion)

- Acute emergency situation
- serious coagulation disorders
- Severe impairment of lungs or cardiovascular system

## 3 Storage

Temperature: -20 °C - +50°C

Relative humidity: 0 -75%, not condensing

Atmospheric pressure: 500 -1600 mbar

Products must be protected from direct sunlight in storage. The storage of the products in their original packaging until their first use is recommended.

## 4 Cleaning and sterilization

Products are supplied in an unsterile state and must be cleaned and sterilised before use.

The electrode handles are sterile packed and must not be cleaned and sterilized before use. They must be disposed of after single use.

Do not clean HF cables, adapters or electrode handles in an ultrasonic bath.

Do not sterilise HF cables, adapters or electrode handles and electrodes in hot air sterilisers.

If a Washing-Desinfection unit is used for cleaning: max 95 °C. The HF cable or handle also may be immersed in a combined cleaning and disinfection solution (e.g. neodisher® MediClean forte, Dr. Weigert Co.). In such a case, the service life may then be reduced due to oxidation of the contacts and crystallisation in the plugs. Please observe the chemical manufacturer's instructions. Then rinse cable, connectors, and handle or electrode thoroughly with demineralized water. Steam-sterilize products at 121°C / 20 min minimum dwell time or 132 resp. 134°C / 4min minimum dwell time (max. 10 min).

## 5 Instructions for Use

Thoroughly examine all products for visible irregularities. Do not use cables with a brittle or defective insulation. In case of doubt check operation by means of a suitable testing device.

### 5.1 Connection to HF unit in monopolar mode

Connect the the foot switch and the "active" HF cable to the HF unit. In order to avoid interference with camera system, do not run HF cables directly in line with camera cables. Never place HF cables directly on the patient's skin. Never lay HF cables in loops. Next, connect the instrument or handle with coagulation electrode, and switch on the HF unit. The electrode with the proper shaft diameter (4 mm or 2.4 mm) must be firmly attached to the handle tip.

### 5.2 Connection to HF unit in bipolar mode

In bipolar HF surgery the coagulation is restricted to the area between the two electrodes of the bipolar instrument. Follow the recommendations regarding patient preparation in monopolar operating mode and the corresponding instructions for the HF unit. When using the bipolar "Auto start" mode, the cable used should not be longer than 4.5 m. Alternatively, use a twin cable with max. length of 5m.

## 6 Warnings

The products may only be used in combination with HEBU accessories and only by clinically instructed and qualified personnel. In particular, the following must be checked: correct attachment of all HF cables and plugs, correct activation of current by the finger - or footswitches, undamaged insulation of the HF cables and instrument, and check that the distal end of the electrode is clean and undamaged. Before operating, check that the signal transmission of the monitors works without interference (e.g. without noise); this check is made by pressing the yellow button (CUT) and the blue button (COAG).

Always unplug the HF-cable by holding on the connector. **Never pull the cable itself**, as this may lead to defects. Cables with defective/brittle insulation or with a broken electrical wire (due to severe bending) might cause a fire. Active HF instruments and HF handles must not be placed on the patient in order to avoid injuries in case of involuntary HF activation or activated "Auto Start" mode.

**Endogenous burns** are those caused by high current density in the patient's tissue. Possible reasons are: The patient is inadvertently positioned in contact with electrically conductive parts. Direct contact between skin areas and HF cables may lead to electrical capacitance which in turn may cause burns.

**Exogenous burns** are those caused by heat of ignited fluids or gases. They may also be caused by explosions. Possible reasons include: Ignition of skin cleaning agents and disinfectants, ignition or narcotic gases etc.

**Cardiac pacemakers** may be damaged by electrosurgical current. Consult a cardiologist before operating. Never use electrosurgical current on outpatients fitted with pacemakers.

## 7 Safety Remarks

	Check the isolation on the electrode and the shaft extensions for damage before each use.
	If damaged, do not use these products!
	Do not use in the presence of combustible or explosive materials!
	Never lay instruments on the patient or in his/her direct vicinity!
	HEBUmedical expressly warns against modifying the product. Any modification exempts HEBUmedical in from any and all liability. .
	Protect these products from any form of mechanical damage! Do not throw! Do not: use force!
	Defective HF cables or electrode handles may not be repaired. Always replace a defective part with a new one.

## 8 Servicing and manufacturer address

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**Our vast range of products covers over 10 000 various instruments.  
They can be found in Internet or request our catalog.**

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**HEBUmedical GmbH**