HF- Monopolare Elektroden

HF monopolar electrodes

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Seite 3

Operating Manual

Page 8

DEUTSCH

ENGLISH

Inhalt

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Contents

1	Used Symbols	9
2	Scope	9
2.1	Inteded Use	9
2.2	Contraindication	10
3	Storage	10
4	Cleaning and sterilization	10
5	Instructions for Use	10
5.1	Connection to HF unit in monopolar mode.....	10
6	Warnings	11
7	Safety Remarks	12
8	Servicing and manufacturer address	12

1 Used Symbols

Symbol	Definition
	CE-marking
	Attention
	Keep away from sunlight
	Store at dry place
	No intervention on patients with cardiac pacemakers.
	Manufacturer

2 Scope

The use of HF monopolar electrodes takes place in HF surgery for coagulation and tissue cutting. Among others:

- transection of lips
- Opening of acute abscesses
- Periodontal disease
- For the removal of tumors and fistulas
- To expose retained teeth

2.1 Inteded Use

HF electrodes, monopolar; HF argon monopolar electrodes; HF monopolar instrument with exchangeable electrodes: A monopolar electrical conductor (blade or tip) connected to a mating electrical handpiece for establishing an active electrical connection between the output terminals of an electrosurgical generator and a patient. The delivered energy is used to perform lesions (cutting and / or cauterizing) at a desired location on the patient. It is commonly used in various medical disciplines. It is a reusable product.

HF coagulation suction cannula: A monopolar suction tube for aspirating fluids and coagulating soft tissue during surgical procedures. It is a reusable product.

Maximum electrical capacity: 4,3kV_p.

HF coagulation suction tube

A Monopolar suction tube for suction of fluids and for coagulating soft tissue during surgery.

Maximum electrical capacity: 4,3kV_p

When used as intended, these products are more or less subject to wear, depending on the intensity of use. This kind of wear is caused by technical factors and cannot be avoided. Therefore carry out a thorough visual examination and a function test prior to each application of these products. If a new product has externally visible defects or does not operate according to the instructions, contact the manufacturer or his responsible representative immediately.

2.2 Contraindication

The HF monopolar electrodes should not be used in conjunction with the central nervous or circulatory system. Tube ligation is not a safe application. Patients with cardiac pacemakers and in-vivo cardiac defibrillators have limitations in their use. RF Monopolar Electrodes should not be used in the immediate vicinity of bone, tooth root, and metal.

3 Storage

Temperature: -20 °C - +50°C

Relative humidity: 0 -75%, not condensing

Atmospheric pressure: 500 -1600 mbar

Products must be protected from direct sunlight in storage. The storage of the products in their original packaging until their first use is recommended.

4 Cleaning and sterilization

Products are supplied in an unsterile state and must be cleaned and sterilised before use.

Do not clean HF cables, adapters or electrode handles in an ultrasonic bath.

Do not sterilise HF cables, adapters or electrode handles and electrodes in hot air sterilisers.

If a Washing-Desinfection unit is used for cleaning: max 95 °C. The HF cable or handle also may be immersed in a combined cleaning and disinfection solution (e.g. neodisher® MediClean forte, Dr. Weigert Co.). In such a case, the service life may then be reduced due to oxidation of the contacts and crystallisation in the plugs. Please observe the chemical manufacturer's instructions. Then rinse cable, connectors, and handle or electrode thoroughly with demineralized water. Steam-sterilize products at 121°C / 20 min minimum dwell time or 132 resp. 134°C / 4min minimum dwell time (max. 10 min).

5 Instructions for Use

Thoroughly examine all products for visible irregularities. Do not use cables with a brittle or defective insulation. In case of doubt check operation by means of a suitable testing device.

5.1 Connection to HF unit in monopolar mode

First fix the neutral electrode ("return plate") to the patient, preferably on the upper arm or on the thigh. Ensure that the skin is free from hair and grease. If a reusable neutral electrode (return plate) is used, do not apply conductive gel. The current paths inside the patient's body should be short and must run diagonally. Current paths must never run in a transverse direction through the body or across the thorax. The patient must be insulated against all electrically conductive parts. Ground the operating table. Place the patient on a dry electrically insulating layer. Prevent skin areas touching

each other. Place dry gauze in between. Then connect the neutral electrode, the foot switch and the “active” HF cable to the HF unit. In order to avoid interference with camera system, do not run HF cables directly in line with camera cables. Never place HF cables directly on the patient’s skin. Never lay HF cables in loops. Next, connect the instrument or handle with coagulation electrode, and switch on the HF unit. The electrode with the proper shaft diameter (4 mm or 2.4 mm) must be firmly attached to the handle tip. For handles with finger switches:

Yellow switch = actuation of cutting current

Blue switch = actuation of tissue coagulation

Always refer to the instruction manual of the HF unit and to the general information it provides relating to electrosurgical procedures.

6 Warnings

The products may only be used in combination with HEBU accessories and only by clinically instructed and qualified personnel. Electrodes must not be twisted, as this may cause damage to the material and the electrode might break. Power output of the HF unit: Always select the lowest possible output of the HF unit for any procedure. If the coagulation capacity of the electrode is below normal, do not increase the output of the HF unit without a thorough preliminary check. In particular, the following must be checked: correct attachment of all HF cables and plugs, correct activation of current by the finger - or footswitches, undamaged insulation of the HF cables and instrument, and check that the distal end of the electrode is dean and undamaged. Before operating, check that the signal transmission of the monitors works without interference (e.g. without noise); this check is made by pressing the yellow button (CUT) and the blue button (COAG).

Always unplug the HF-cable by holding on the connector. **Never pull the cable itself**, as this may lead to defects. Cables with defective/brittle insulation or with a broken electrical wire (due to severe bending) might cause a fire. Active HF instruments and HF handles must not be placed on the patient in order to avoid injuries in case of involuntary HF activation or activated “Auto Start” mode.

For cleaning the active electrode, you must deactivate the “Auto Start” mode (only in bipolar applications) or disconnect the instrument from the generator. Directly after activation the tips could be hot. This may lead to burns.

Endogenous burns are those caused by high current density in the patient’s tissue. Possible reasons are: The patient is inadvertently positioned in contact with electrically conductive parts. Direct contact between skin areas and HF cables may lead to electrical capacitance which in turn may cause burns.

Exogenous burns are those caused by heat of ignited fluids or gases. They may also be caused by explosions. Possible reasons include: Ignition of skin cleaning agents and disinfectants, ignition or narcotic gases etc.

Cardiac pacemakers may be damaged by electrosurgical current. Consult a cardiologist before operating. Never use electrosurgical current on outpatients fitted with pacemakers.

7 Safety Remarks


	Check the isolation on the electrode and the shaft extensions for damage before each use.
	Do not rebend needles, loops or knives. Spatulas can be bent.
	If damaged, do not use these products!
	Do not use in the presence of combustible or explosive materials!
	Never lay instruments on the patient or in his/her direct vicinity!
	HEBUmedical expressly warns against modifying the product. Any modification exempts HEBUmedical in from any and all liability. .
	Protect these products from any form of mechanical damage! Do not throw! Do not: use force!
	Defective HF cables or electrode handles may not be repaired. Always replace a defective part with a new one.

8 Servicing and manufacturer address

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**Our vast range of products covers over 10 000 various instruments.
They can be found in Internet or request our catalog.**

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HEBUmedical GmbH