Diagnoseinstrumente / Diagnostic Instruments

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Seite 3

Operating Manual

Page 11

DEUTSCH

ENGLISH

Inhalt

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Contents

1	Symbols	11
2	Introduction	11
3	Scope	11
3.1	Inteded Use	11
3.2	Contraindications.....	12
4	Warning.....	12
5	Handling	12
6	Preparation.....	13
6.1	Reutilization restrictions	13
6.2	Information on instrument preparation	13
6.3	Preparation at the place of use	13
6.4	Ultrasound bath (optional)	13
6.5	Manual cleaning	13
6.6	Mechanical cleaning.....	14
6.7	Drying.....	15
7	Sterilization.....	15
8	Maintenance, inspection	16
9	Packaging	16
10	Storage.....	16
11	Warranty / Repair	16
12	Manufacturer and servicing address.....	16

1 Symbols

Symbol	Definition
	CE marking
	Attention
	Manufacturer

2 Introduction

By purchasing this instrument, you are now the owner of a high-quality product whose use and correct handling are described in the following. In order to minimize possible risks to patients and users, please observe these instructions carefully. Use, disinfection, cleaning and sterilization may only be performed by suitably trained specialist personnel.

3 Scope

The use of measuring Instruments is mainly in the following areas:

- Surgery
- Orthopedics
- ENT

3.1 Inteded Use

Finger-Goniometer; Goniometer(Flexometer); Testing and measuring instrument: An instrument for comparative measurements of the degree of rotation, e.g. from the eyeball or shaft of a long bone. The instrument is not calibrated and can be reused.

Dynamometer: A product to investigate and determine the level of vibration needed to determine the vibratory threshold of perception on the skin anywhere on the body. It is used for neuropathological diagnoses, early detection of dysfunction, e.g. in the case of suspected diabetes, patients with neurotoxic mixtures, unfavorable working conditions or car accidents. It is a reusable product.

Neurological Pinwheel: A manual sensiometer for determining the tactile sensation of a patient by distinguishing the distance of rods mounted on a disc. It typically consists of a disc with protruding rods that differ in diameter. The disc is mounted on a rotating pin. The disc is manually moved on the skin of the patient via an integrated handle. The rods can be pointed to test the pain sensation. It is a reusable product.

Pelvimeters: A medical device for comparison measurement of the pool size. These are usually the inner / outer dimensions of the pelvis, its diameter and its capacity. It is commonly used in gynaecological procedures / procedures. The instrument is not calibrated and can be reused.

Laser-Instruments; Width ruler; Depth ruler: A surgical instrument used to determine the depth of a cavity (usually a drilled hole). It usually consists of a movable (sliding) central piece (a pin) whose point is inserted into the cavity; if the tip touches the bottom of the cavity, the value can be read visually on a marked scale on the side of the instrument. It is a comparative measuring instrument and is often used to determine the depth of drill holes in bones so that the surgeon can determine the length of bone screws to use thereafter. It is made of corrosion-resistant materials (eg: stainless steel, titanium). It is a reusable instrument.

Marker graduated: An ophthalmic instrument consisting of two legs pivotally anchored at a point used to measure the diameter, length, angle and thickness of the eye. The instrument is not calibrated and is used for comparative measurements. It is a reusable product.

3.2 Contraindications

There are no known contraindications.

4 Warning

	Medical products are delivered in a non-sterile condition and must be cleaned, disinfected and sterilized prior to their initial use.
	The use of faulty instruments is in principle forbidden and they have to go through the whole cleaning process before return.
	Please take into consideration that through higher power a bigger damage of the tissue can result: f.e. on forceps: the power at the end of the jaw is higher than at the tip of the jaw
	Please observe the additional information enclosed with the products.
	Remove all protective sleeves and films prior to first using or preparation for use.
	The safe combination of different products or of products with implants must be reviewed prior to clinical application by the user.
	Avoid improper throwing or dropping of instruments
	To prevent all contact corrosion, instruments with damaged surfaces must be separated immediately.
	If used on patients with Creutzfeldt-Jakob disease or HIV infection, we are unable to accept any responsibility for reuse.

5 Handling

The type of treatment must be determined in each individual case by the surgeon in cooperation with the internist and the anaesthetist.

For operational use in various surgical disciplines must be done by appropriately trained and qualified personnel.

6 Preparation

The person in charge of preparatory treatment is responsible for ensuring that the treatment is duly carried out using the relevant equipment, materials and personnel in the treatment facility and so achieves the desired result. This necessitates validation and routine monitoring of the process used. We urge you to take note of the national regulations dealing with instrument preparation.

6.1 Reutilization restrictions

Frequent repeat preparatory treatment has minimal effects on the product. The end of the product life is normally determined by wear and damage due to use

6.2 Information on instrument preparation

- Use cleaning and/or disinfection agents with a pH-value within 9-10. Please observe manufacturer instructions regarding dosage, exposure time and renewal of solutions.
- Do not use hard brushes or coarse abrasive cleaners.
- Never leave instruments in cleaning or disinfection agents for longer than the specified time.
- Only used demineralized water for rinsing.
- Rinse and dry carefully through channels and pipes.
- Sensitive instruments must be cleaned in a storage or clamping fixture.
- Observe manufacturer instructions of cleaning – and sterilizing equipment.

6.3 Preparation at the place of use

Directly after using remove coarse dirt of the instruments and rinse out the working cannulas. Do not use fixing agents or hot water ($> 40^{\circ}\text{C}$), as this results in residues becoming fixed and can affect the success of the subsequent cleaning operation

Dismantle and/or open instruments as far as possible. Within short time after use the instruments clean the instruments for reducing a drying of the residues. **Use cleaning and/or disinfection agents with a neutral pH-value < 10 .**

This enables an easier cleaning. If instruments come into contact with corroding medicines or cleaning agents, wash these up with water immediately after use.

6.4 Ultrasound bath (optional)

All instruments must be opened, dismantled and any cavities rinsed through.

Place instruments in the screen basket in such a way that overlaps and contact between instruments are avoided. Add cleaning

agent to the water and adjust the temperature of the solution in line with the cleaning agent manufacturer's instructions.

The cleaning in the ultrasound bath should be at 35-40 kHz, 5 minutes at least.

Subsequently rinse instruments including all cavities before cleaning and disinfection.

Medicine products which are possessing a bad ultrasound transmission, e.g. soft material are not usable for the ultrasound bath.

6.5 Manual cleaning

We recommend using lint-free soft cloths, plastic brushes and/or cleaning guns for manual cleaning. The cleaning agent should not have any fixing properties. The cleaning fluid must be changed daily.

Hints:

- If possible, apply the cleaning agent into cavities at pressure and carefully brush out.
- Be careful with the cleaning of barriers, slits, tight entries and joints etc.
- It has to be secured that instruments with joints will be cleaned in an opened and closed position.

- Use adequate brushes for the cleaning of cannulas, working channels and sack holes
- Flexible (complex) instruments with unvisible surface and chuck have to be cleaned manually.
- We recommend for all push shafts, -and pipe shafts instruments and instruments whose surfaces are on top of each others during the cleaning (f.e. bone forceps and gouge forceps) a manually cleaning for an optimal cleaning result without residues

After that rinse the instruments with fully demineralized water extensively. Subsequently dip the instrument in disinfectant and use suitable cleaning instruments such as synthetic soft brushes. Observe the cleaning agent manufacturer's instructions. At the end of the cleaning process, repeat the rinsing action with fully demineralized water and dry the instruments completely.

6.6 Mechanical cleaning



On the basis of international standards (EN ISO 15883) and national directives, only validated machine cleaning and disinfection methods may be used. For the mechanical cleaning we recommend a standard programme for surgical instruments, f.e. instruments from Miele.

Observe the following by loading:

- Place the dismantled/opened instruments securely in the tray.
- Instruments with openings and gaps have to be faced down with the opened side so that they can be cleaned and no water of the cleaning process is collecting inside them. If available use balanced devices for rinsing

- Place the instruments with joints in an opened position into the cleaning, -and disinfection machine.
- If needed use an adapter for the cleaning
- Do not overload trays, avoid creating any overlaps.

Preliminary rinsing (cold, if applicable fully demineralized water without additives) is followed by chemical.

The chemical cleaning should take place at 40°C -60°C for at least 5 minutes.

We recommend products with a **pH-value within 9-10**, e.g. MediClean forte from Dr.Weigert. The cleaning agents used should be selected depending on the material and properties of the instruments and in accordance to national regulations: If there is a high chloride concentration in the water, pitting and tension crack corrosion can occur on the instruments. The occurrence of this type of corrosion is minimized by using alkaline cleaning agents and demineralized water. By adding an acid-based neutralization medium, the rinsing off of alkaline cleaning agent residues is facilitated during the first intermediate rinsing process (warm or cold water).

In order to prevent the formation of deposits, it is advisable to use neutral cleaners where the water quality is unfavourable. After the second intermediate rinsing process, thermal disinfection takes place.

The thermal disinfection should take place at temperatures of between 80 and 95°C, with an exposure time as outlined in EN ISO 15883.

After the finishing of the programme take the good out of the machine because corrosion can arise if the instrument remains in the machine.

6.7 Drying

Ensure adequate drying by the cleaning and disinfection device or using other suitable measures.

The drying temperature should not exceed 95°C to avoid material-related ageing processes.

7 Sterilization

Prior to sterilization, products must undergo cleaning and disinfection, be rinsed off without residue using demineralized water and subsequently dried. HEBUmedical recommends using a validated steam sterilization process (e.g. sterilizer in compliance with EN 285 and validated in accordance with DIN EN ISO 17665-1).

On using the fractionated vacuum method, **sterilization** must be performed with at least **134°C (USA 132° C)** with a **minimum dwell period of 4 minutes** and subsequent vacuum drying.

The vapour must be free of ingredients, recommended limiting values of feed water and vapour condensate are determined through EN 285.

Other sterilization processes are compatible but not validated from HEBUmedical.

When loading, observe the recommended total weight. After the sterilization, check the sterile product packaging for damage, and inspect the sterilization indicators

8 Maintenance, inspection

After cooling to room temperature, the instruments must be visually inspected for protein residues and other contamination, paying particular attention to cavities, blocks, inclusions, pipes, and other inaccessible areas. Instruments which are not free of residues must be returned for a complete retreatment process.

Carry out the functional check mentioned above. Instruments with stains, which are blunt, bent, no longer function or which are otherwise damaged must be segregated!

9 Packaging

Compliant packaging of products for sterilization in line with ISO 11607. Packaging used must be suitable for the instruments and protect them from microbiological contamination during storage. The seal must not be under tension. HEBUmedical recommends container for sterilized packaging.

10 Storage

Store products in a dry, clean and dust-free environment at moderate temperatures from 5°C to 40°C. Protect from the effects of the sun's rays and artificial light. Warranty / Repair

11 Warranty / Repair

Our products are manufactured from high-grade materials and carefully checked prior to dispatch. However, even if used properly in accordance with their intended purpose they are subject to a greater or lesser degree of wear depending on their intensity of use.

This wear is technically induced and unavoidable.

Should faults occur independently of wear, please contact our customer services. Defective products should no longer be used.

They must undergo the complete preparatory treatment process before being returned.

12 Manufacturer and servicing address

HEBUmedical GmbH
Badstraße 8
78532 Tuttlingen / Germany
Tel. +49 7461 94 71 - 0
Fax +49 7461 94 71 - 22
eMail: service@HEBUmedical.de
Web: www.HEBUmedical.de

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**Our vast range of products covers over 10 000 various instruments.
They can be found in Internet or request our catalog.**

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