Einfache Schneidinstrumente / Simple cutting instruments

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Operating Manual
Page 14

DEUTSCH

ENGLISH

Inhalt

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Contents

| | | |
|-----|---|----|
| 1 | Symbols | 14 |
| 2 | Introduction | 14 |
| 3 | Scope..... | 14 |
| 3.1 | Inteded Use | 14 |
| 3.2 | Contraindications..... | 16 |
| 4 | Warning..... | 17 |
| 5 | Handling..... | 17 |
| 6 | Preparation | 17 |
| 6.1 | Reutilization restrictions..... | 17 |
| 6.2 | Information on instrument preparation | 18 |
| 6.3 | Preparation at the place of use | 18 |
| 6.4 | Ultrasound bath (optional) | 18 |
| 6.5 | Manual cleaning | 18 |
| 6.6 | Mechanical cleaning | 19 |
| 6.7 | Drying..... | 20 |
| 7 | Sterilization | 20 |
| 8 | Maintenance, inspection..... | 20 |
| 9 | Packaging | 21 |
| 10 | Storage | 21 |
| 11 | Warranty / Repair | 21 |
| 12 | Manufacturer and servicing address..... | 21 |

1 Symbols

| Symbol | Definition |
|---|--------------|
|  | CE marking |
|  | Attention |
|  | Manufacturer |

2 Introduction

By purchasing this instrument, you are now the owner of a high-quality product whose use and correct handling are described in the following. In order to minimize possible risks to patients and users, please observe these instructions carefully. Use, disinfection, cleaning and sterilization may only be performed by suitably trained specialist personnel.

3 Scope

The use of simple cutting instruments is mainly in the following areas:

- abdominal surgery
- Gynecology
- pathology
- Urology

3.1 Inteded Use

Knives

Amputating Knives: A heavy-duty surgical instrument with a handle and a single-sided or double-sided cutting blade of various shapes and sizes for use in limb amputations. It is a reusable instrument.

Flap knife; Phalangeal knife; Interosseous knife; Post mortem knives; Resection knives; Autopsy knife; Micro Vessel Knives; Micro-laryngeal Instruments; Knives; Mucosa Knife; Knives, others: A surgical instrument with a handle and a scalpel blade (not interchangeable) used by doctors to cut or prepare tissues with downward pressure. The instrument is made of stainless steel and is reusable. For a scalpel with replaceable blade, see the corresponding generic product group.

Plaster-of-Paris Knives: A surgical hand instrument with a short, flat, strong cutting blade with a sharp edge for cutting or shortening plaster casts. It is usually a one-piece product made of stainless steel. It is a reusable product. Rhinoplastic and Nasal Knives; Tonsil Knives, Dissectors, Enucleators, Retractors; Tympanum Knives; Trigeminal Knives; A surgical hand instrument for cutting anatomical conditions during a general ENT procedure. It is made of stainless steel and has a handle in various designs, which merges into a narrow working end with a narrow, sharp blade. It is a reusable instrument.

Antrum Punch; Sphenoid Punches: A surgical instrument for the removal of tissue or bone samples during ear, nose and throat (ENT) procedures. The instrument can be equipped with a ring handle or pistol grip. The ring handles or pistol grips have extended shafts provided with punch-like or mechanical locking mechanisms at the distal end. The instrument is a reusable instrument.

Adenoid Curettes: A surgical hand instrument for cutting and excision of lymphoid tissue during an adenoidectomy. It is usually performed as a long, narrow instrument with a handle and a concave, spoon-like working part. The working part has a sharp edge for removing the tissue without injuring the surrounding muscles. It is usually made of stainless steel. It is a reusable instrument.

Swivel Knives: A surgical hand instrument for cutting and excision of lymphoid tissue during an adenoidectomy. It is usually performed as a long, narrow instrument with a handle and a concave, spoon-like working part. The working part has a sharp edge for removing the tissue without injuring the surrounding muscles. It is usually made of stainless steel. It is a reusable instrument.

Eye Knives; Keratomes: An ophthalmic surgical instrument for placing precise cuts in ocular tissue during ophthalmic surgery on the eye and surrounding structures. It is usually designed as a one-piece instrument with a sharp, single-edged blade at the distal end and a handle at the proximal end. It can have different blade shapes, e.g. phaco (dull or sharp), crescent-shaped, straight, and micro-phaco (for cataract operations). It is usually made of stainless steel and some models may be equipped with diamond blades. It is a reusable instrument.

Knives, others: A surgical hand instrument, usually softly rounded and bent at the working end, to separate the outermost membrane surrounding the brain and spinal cord from the inner membranes. It usually consists of a handle that runs off to a shaft and has a rounded and bent tip at the distal end. It can also be double-ended. It is made of stainless steel and is available in different shapes and sizes. It is a reusable instrument.

Mirco Knives: A surgical instrument for cutting tissue that consists of a long slender metal handle that tapers towards the distal end. The internal radius of the hook has a cutting edge. It is a reusable instrument.

Micro-Instruments: A shaft-like surgical instrument whose distal end has the form of a hook or more parallel hooks for fixing or pulling uterine fibroids during a gynaecological procedure. It is typically a one-piece stainless-steel instrument with a handle in various designs at the proximal end. It is a reusable instrument.

Scalpels

Operating Knives: A surgical instrument with a handle and a scalpel blade (not interchangeable) used by doctors to cut or prepare tissues with downward pressure. The instrument is made of stainless steel and is reusable.

Scalpel Blades: An interchangeable component of a scalpel that works as a cutting edge and is inserted into a suitable grip. Together, it provides a surgical instrument for cutting or dissecting tissues. For scalpel handles, there is a matching generic product group. The blades are delivered sterile and are disposable.

Scalpel Handles: An interchangeable component of a scalpel that works as a handle and picks up a matching blade. Together, it provides a surgical instrument for cutting or dissecting tissues. The instrument can be made of stainless steel or brass. For one-piece scalpels, there is a suitable generic product group. It is a reusable instrument.

Chisel

Bone Chisels and Gouges; Nasal Osteotomes; flat chisel; Rhachiotom: A surgical, chisel-like hand instrument for cutting and / or shaping bones by bumping during orthopedic surgery. It is held by the surgeon who applies a manual force to the proximal end of the instrument via a surgical hammer. The distal end (the cutting or sharp edge) is sharp, often flat, sometimes curved (concave), and usually folded on both sides. It is made of stainless steel and is a reusable instrument.

Sternum chisel: A single-sided tapered blade surgical instrument for cutting and shaping bone during orthopedic surgery. The surgeon holds it in one hand while hitting the proximal end with a surgical gavel in the other hand. The distal end (the cutting edge) is sharp, flat or curved (concave), also called a hollow blade). It is made of high quality stainless steel and may consist of one piece, or have a synthetic handle (such as Tufnol). The instrument is reusable.

Septum Gouges: A surgical instrument with a single sided beveled blade for cutting and shaping nasal bones during ear, nose and throat (ENT) operations. The surgeon holds it in one hand while beating the proximal end of the instrument with a surgical gavel in the other hand. The distal end (the cutting edge) is sharp, but often flat but may be curved (concave). It is typically one piece and made of high quality stainless steel. The instrument is reusable. Others

Dermatom: A surgical hand instrument for cutting thin skin slices for transplantation or excision of small skin lesions. This requires a special cutting blade, which is used for this purpose in the instrument. It is a reusable instrument.

Separators: A surgical hand instrument, usually in spoon shape or rounded at the working end. It is used to separate soft tissue or body structures in interventions of general or plastic surgery. It usually has a handle that continues into a shaft that has a tip at the distal end. The tip may be pointed or flat, sharp or dull, angled or straight on the shaft. The instrument is made of high quality stainless steel and available in various shapes and dimensions. The instrument is reusable.

Enucleators: A surgical hand instrument for the dissection of the tonsils during an ENT procedure. It is usually spoon-like or rounded at the working ends. It is manufactured in various designs, e.g. with a shaft, which merges at the distal end in a rounded, bent or angled tip, or even double-ended. The working end may be tapered, flat, sharp or dull, straight or angled, serrated or smooth. It is made of stainless steel and is available in different sizes and shapes. It is a reusable instrument.

Circumcision Instruments: A surgical instrument for the controlled removal of the penis foreskin during circumcision. It is typically made of high quality stainless steel and may have various shapes, e.g. the bell clamp (bell-shaped clamp). The foreskin is pushed over the bell-shaped protective body and the glans are positioned therein. A screw mechanism is tightened so that the foreskin is compressed in an annular manner and fitted with e.g. a scalpel along the compression ring can be cut off. The instrument is reusable.

Adenotomes: A surgical instrument designed to excise hypertrophic lymphoid tissue in the nasopharynx (i.e., pharyngeal tonsils or polyps) during an adenoidectomy. It is usually a long, slender, manual hand instrument with cutting edges at the working end and a handle with a mechanism for actuating the blades. It is available in various designs and sizes and can be a one-piece instrument or require the insertion of blades into the handle. This is a reusable instrument.

Meniscotomy Knives: A cutting, surgical instrument whose handle is T-shaped to be able to apply more force. The handle merges into a chisel-shaped edge. The instrument is used to cut crescent-shaped parts of the capsule or cartilage area.

Tendon Stripper: A surgical instrument consisting of a handle and a semi-circular, domed working end that ends in a sharp cutting edge. Tendon strippers are used to cut out a piece of ligament, tendon or fascia to transplant them.

3.2 Contraindications

There are no known contraindications.

4 Warning

| | |
|---|--|
|  | Medical products are delivered in a non-sterile condition and must be cleaned, disinfected and sterilized prior to their initial use. |
|  | The use of faulty instruments is in principle forbidden and they have to go through the whole cleaning process before return. |
|  | Please take into consideration that through higher power a bigger damage of the tissue can result: f.e. on forceps: the power at the end of the jaw is higher than at the tip of the jaw |
|  | Please observe the additional information enclosed with the products. |
|  | Remove all protective sleeves and films prior to first using or preparation for use. |
|  | The safe combination of different products or of products with implants must be reviewed prior to clinical application by the user. |
|  | Avoid improper throwing or dropping of instruments |
|  | To prevent all contact corrosion, instruments with damaged surfaces must be separated immediately. |
|  | If used on patients with Creutzfeldt-Jakob disease or HIV infection, we are unable to accept any responsibility for reuse. |

5 Handling

The type of treatment must be determined in each individual case by the surgeon in cooperation with the internist and the anaesthetist.

For operational use in various surgical disciplines must be done by appropriately trained and qualified personnel.

6 Preparation

The person in charge of preparatory treatment is responsible for ensuring that the treatment is duly carried out using the relevant equipment, materials and personnel in the treatment facility and so achieves the desired result. This necessitates validation and routine monitoring of the process used. We urge you to take note of the national regulations dealing with instrument preparation.

6.1 Reutilization restrictions

Frequent repeat preparatory treatment has minimal effects on the product. The end of the product life is normally determined by wear and damage due to use

6.2 Information on instrument preparation

- Use cleaning and/or disinfection agents with a pH-value within 9-10. Please observe manufacturer instructions regarding dosage, exposure time and renewal of solutions.
- Do not use hard brushes or coarse abrasive cleaners.
- Never leave instruments in cleaning or disinfection agents for longer than the specified time.
- Only used demineralized water for rinsing.
- Rinse and dry carefully through channels and pipes.
- Sensitive instruments must be cleaned in a storage or clamping fixture.
- Observe manufacturer instructions of cleaning – and sterilizing equipment.

6.3 Preparation at the place of use

Directly after using remove coarse dirt of the instruments and rinse out the working cannulas. Do not use fixing agents or hot water ($> 40^{\circ}\text{C}$), as this results in residues becoming fixed and can affect the success of the subsequent cleaning operation

Dismantle and/or open instruments as far as possible. Within short time after use the instruments clean the instruments for reducing a drying of the residues. **Use cleaning and/or disinfection agents with a neutral pH-value <10 .**

This enables an easier cleaning. If instruments come into contact with corroding medicines or cleaning agents, wash these up with water immediately after use.

6.4 Ultrasound bath (optional)

All instruments must be opened, dismantled and any cavities rinsed through.

Place instruments in the screen basket in such a way that overlaps and contact between instruments are avoided. Add cleaning

agent to the water and adjust the temperature of the solution in line with the cleaning agent manufacturer's instructions.

The cleaning in the ultrasound bath should be at 35-40 kHz, 5 minutes at least.

Subsequently rinse instruments including all cavities before cleaning and disinfection.

Medicine products which are possessing a bad ultrasound transmission, e.g. soft material are not usable for the ultrasound bath.

6.5 Manual cleaning

We recommend using lint-free soft cloths, plastic brushes and/or cleaning guns for manual cleaning. The cleaning agent should not have any fixing properties. The cleaning fluid must be changed daily.

Hints:

- If possible, apply the cleaning agent into cavities at pressure and carefully brush out.
- Be careful with the cleaning of barriers, slits, tight entries and joints etc.
- It has to be secured that instruments with joints will be cleaned in an opened and closed position.
- Use adequate brushes for the cleaning of cannulas, working channels and sack holes
- Flexible (complexe) instruments with unvisible surface and chuck have to be cleaned manually.
- We recommend for all push shafts, -and pipe shafts instruments and instruments whose surfaces are on top of each others during the cleaning (f.e.bone forceps and gouge forceps) a manually cleaning for an optimal cleaning result without residues

After that rinse the instruments with fully demineralized water extensively. Subsequently dip the instrument in disinfectant and use suitable cleaning instruments such as synthetic soft brushes.



Observe the cleaning agent manufacturer's instructions. At the end of the cleaning process, repeat the rinsing action with fully demineralized water and dry the instruments completely.

6.6 Mechanical cleaning



On the basis of international standards (EN ISO 15883) and national directives, only validated machine cleaning and disinfection methods may be used. For the mechanical cleaning we recommend a standard programme for surgical instruments, f.e. instruments from Miele.

Observe the following by loading:

- Place the dismantled/opened instruments securely in the tray.
- Instruments with openings and gaps have to be faced down with the opened side so that they can be cleaned and no water of the cleaning process is collecting inside them. If available use balanced devices for rinsing.

- Do not overload trays, avoid creating any overlaps.

Preliminary rinsing (cold, if applicable fully demineralized water without additives) is followed by chemical.

The chemical cleaning should take place at 40°C -60°C for at least 5 minutes.

We recommend products with a **pH-value within 9-10**, e.g. MediClean forte from Dr.Weigert. The cleaning agents used should be selected depending on the material and properties of the instruments and in accordance to national regulations: If there is a high chloride concentration in the water, pitting and tension crack corrosion can occur on the instruments. The occurrence of this type of corrosion is minimized by using alkaline cleaning agents and demineralized water. By adding an acid-based neutralization medium, the rinsing off of alkaline cleaning agent residues is facilitated during the first intermediate rinsing process (warm or cold water).

In order to prevent the formation of deposits, it is advisable to use neutral cleaners where the water quality is unfavourable. After the second intermediate rinsing process, thermal disinfection takes place.

 The thermal disinfection should take place at temperatures of between 80 and 95°C, with an exposure time as outlined in EN ISO 15883.

After the finishing of the programme take the good out of the machine because corrosion can arise if the instrument remains in the machine.

6.7 Drying

Ensure adequate drying by the cleaning and disinfection device or using other suitable measures.

 The drying temperature should not exceed 95°C to avoid material-related ageing processes.

7 Sterilization

Prior to sterilization, products must undergo cleaning and disinfection, be rinsed off without residue using demineralized water and subsequently dried. HEBUmedical recommends using a validated steam sterilization process (e.g. sterilizer in compliance with EN 285 and validated in accordance with DIN EN ISO 17665-1).

 On using the fractionated vacuum method, **sterilization** must be performed with at least **134°C (USA 132° C)** with a **minimum dwell period of 4 minutes** and subsequent vacuum drying.

The vapour must be free of ingredients, recommended limiting values of feed water and vapour condensate are determined through EN 285.

Other sterilization processes are compatible but not validated from HEBUmedical.

When loading, observe the recommended total weight. After the sterilization, check the sterile product packaging for damage, and inspect the sterilization indicators.

8 Maintenance, inspection

After cooling to room temperature, the instruments must be visually inspected for protein residues and other contamination, paying particular attention to cavities, blocks, inclusions, pipes, and other inaccessible areas. Instruments which are not free of residues must be returned for a complete retreatment process.

Carry out the functional check mentioned above. Instruments with stains, which are blunt, bent, no longer function or which are otherwise damaged must be segregated!

9 Packaging

Compliant packaging of products for sterilization in line with ISO 11607. Packaging used must be suitable for the instruments and protect them from microbiological contamination during storage. The seal must not be under tension. HEBUmedical recommends container for sterilized packaging.

10 Storage

Store products in a dry, clean and dust-free environment at moderate temperatures from 5°C to 40°C. Protect from the effects of the sun's rays and artificial light.

11 Warranty / Repair

Our products are manufactured from high-grade materials and carefully checked prior to dispatch. However, even if used properly in accordance with their intended purpose they are subject to a greater or lesser degree of wear depending on their intensity of use.

This wear is technically induced and unavoidable.

Should faults occur independently of wear, please contact our customer services. Defective products should no longer be used.

They must undergo the complete preparatory treatment process before being returned.

12 Manufacturer and servicing address

HEBUmedical GmbH
Badstraße 8
78532 Tuttlingen / Germany
Tel. +49 7461 94 71 - 0
Fax +49 7461 94 71 - 22
eMail: service@HEBUmedical.de
Web: www.HEBUmedical.de

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**Our vast range of products covers over 10 000 various instruments.
They can be found in Internet or request our catalog.**

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