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Tube shaft instruments / Instruments de tuyaux / Instrumentos de tubería / Strumenti a tubo / Тръбни инструменти / Rørinstrumenter / Toruinstrumentid / Putkensoittimet / Όργανα σωλήνων / Instrumenti cijevi / Cauruļu instrumenti / Vamzdžių instrumentai / Pijpinstrumenten / Rørskafinstrumenter / Instrumenty rurowe / Instrumentos de tubo / Instrumente de țevă / Rörinstrument / Potrubné nástroje / Cevni instrumenti / Potrubní nástroje / Cső műszerek

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Seite 5

Operating Manual

Page 23

Mode d'emploi

Page 39

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Página 57

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Página 75

DEUTSCH

ENGLISH

FRANÇAIS

ESPAÑOL

ITALIANO

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Table of contents

1	Symbol descriptions	24
2	Introduction	25
3	Intended use	25
	3.1 Intended purpose	25
	3.2 Indication	27
	3.3 Contraindications	28
	3.4 Patient target group	28
4	Warnings	29
5	Handling	30
6	Preparation	30
	6.1 Information on instrument preparation	30
	6.2 Preparation at the place of use	31
	6.3 Ultrasound bath (optional)	31
	6.4 Manual cleaning	31
	6.5 Mechanical cleaning	31
	6.6 Drying	34
7	Maintenance, inspection	35
	7.1 Functional check	35
8	Sterilization	36
	8.1 Packaging	36
9	Lifetime	36
10	Storage	37
11	Warranty / Repair	37
12	Disposal	37
13	Service and manufacturer address	38

1 Symbol descriptions

Symbols	Definition
	CE-labelling
	Attention
	Validated Parameters
	Manufacturer
	Lot-description
	Reference code
	Medical device / FDA Prescription device
	Medical device
	Non sterile
	Keep away from sunlight
	Dry storage required
	(Electronic) instruction for use

2 Introduction

By purchasing this instrument, you are now the owner of a high-quality product whose use and correct handling are described in the following.

In order to minimize possible risks to patients and users, please observe these instructions carefully. Use, disinfection, cleaning and sterilization may only be performed by suitably trained specialist personnel.

Our products are exclusively intended for professional use by appropriately trained and qualified personnel and may only be acquired by them.

3 Intended use

3.1 Intended purpose

Trocars; Breast Enlarge-ment Sets	A surgical instrument with a sharp pyramidal or conical tip used to pierce body cavity. It is usually assembled and used with a compatible sleeve, thereby allowing insertion of the instrument. After piercing, the trocar is withdrawn and provides a working channel into the body cavity. It is a reusable instrument intended for transient use. It is not intended for use in direct contact with the heart, central circulatory system, or central nervous system.
Haemorrhoidal Ligators; rectal snares	A surgical instrument used to prevent hemorrhoids. A ligator usually consists of a long shaft or cannula with a clamp or handle at the distal end to hold the ligature. The proximal end consists of a handle with moving parts that mechanically actuate the distal clamp or handle. It is a reusable instrument intended for transient use.
Rectal Biopsy Specimen Forceps; Biopsy Punches	A surgical instrument for the general removal of biopsies of tumours and other tissues during surgical intervention for histopathologic investigation. It is typically made of high quality stainless steel and it is a long, inelastic instrument whose distal ends are furnished with two formed cutting shells with sharpened edges, one of them fits exactly in the other one, or scissor-like shells with a raw edge, which offer the cut of the biopsy. These ones are controled due to a scissor-like ring handle on the proximal end. It is inserted over a man-made or a natural body opening in body cavities. The instrument is available in a large number of sizes. It is a reusable instrument intended for transient use. It is not intended for use in direct contact with the heart, central circulatory system, or central nervous system.
Cotton Applicators	A surgical product for use in combination with a suitable endoscope. It is used to bring a piece of absorbent gauze or fabric into the surgical site. It is a reusable instrument intended for transient use.
Bone Shears; Rib Shears	A surgical tool to cut bone or cartilage tissue during orthopedic surgery. It consists of two movable blades, which are usually provided with finger and thumb ring handles. The distal end of the leaves can be designed differently. It is a reusable instrument intended for transient use.

<p>Ligature Guides; Ligature Conductors; Tendon passing instrument; Wire Strut Guide</p>	<p>A surgical instrument for guiding wire or ligatures through tissue. Typical guides may have flexible or rigid working parts with hooks, clips or claws for holding the material while passing through the tissue. The working end can have a button or an eye. It is a reusable instrument intended for transient use. It is not intended for use in direct contact with the heart, central circulatory system, or central nervous system.</p>
<p>Trachea Trocars</p>	<p>A rigid surgical instrument with a sharp pyramidal or conical tip used to puncture body cavities during ear, nose and throat (ENT) surgery. It can be used in conjunction with a compatible sleeve into which it is inserted. Following puncture by the trocar located in the sleeve, the trocar is withdrawn, thereby releasing a working channel through the sleeve into the body cavity. It is a reusable instrument intended for transient use.</p>
<p>Scissors (Tube shaft)</p>	<p>A general surgical instrument for cutting anatomical structures during an ENT procedure at e.g. the bronchi, esophagus, trachea, larynx, pharynx or nose. It has a scissor-type design with ring handles, a tube shaft and a working end with different mouth designs, e.g. straight, angled, or curved. It is made of stainless steel. It is a reusable instrument for transient use.</p>
<p>Goiter forceps; throat forceps; laryngoscopy instruments; laryngeal polyp forceps; nose and ear tweezers; nasal septum forceps; nose pliers; Ear forceps polyp forceps; Sharp nose pliers; septum forceps</p>	<p>A general surgical instrument for transient gripping, holding or manipulating anatomical structures during an surgical ENT procedure at e.g. the bronchi, esophagus, trachea, larynx, pharynx, nose, or ear. It has a scissor-type design with ring handles, a tube shaft and a working end with different mouth designs, e.g. straight, angled, or curved with teeth or ridges to improve retention. It is made of stainless steel. It is a reusable instrument intended for transient use</p>
<p>Needle holders with T/C jaws; Needle holders; Micro needle holders; Needle holders (tube shaft)</p>	<p>A surgical hand-held instrument for gripping surgical needles during a suture to guide the needle and suture through tissue. It is made of stainless steel and has a short-ridged jaw part to allow for improved retention of the needle. It is available in many sizes and designs, e.g. 1) a self-holding, scissor-like design with ring handles; 2) a tweezer-like design with curved, hinged handles that are compressed; and 3) a heavy tong-shaped design. Some models are made with carbide insert in the mouth. It is a reusable instrument intended for transient use.</p>
<p>Prostatic Dilators</p>	<p>A instrument for examining and dilating the urethra. It is a reusable, invasive instrument with respect to body orifices, which are not intended for connection to an active device and for transient use.</p>
<p>Grasping instruments</p>	<p>An instrument for grasping and removing IUDs. It is a reusable, invasive instrument with respect to body orifices, which are not intended for connection to an active device and for transient use.</p>
<p>Sphenoid Punches; Ethmoid Forceps; Antrum Punch</p>	<p>A surgical instrument for the removal of tissue or bone samples during ear, nose and throat (ENT) procedures. The instrument can be equipped with a ring handle or pistol grip. The ring handles or pistol grips have extended shafts provided with punch-like or mechanical locking mechanisms at the distal end. It is a reusable instrument intended for transient use.</p>

Nasal Scissors (Tube shaft); Septum Scissors; Nasal Scissors	A surgical instrument for cutting tissue during an ENT procedure or during plastic surgery on the nose and its adjacent areas. The scissors consists of two movable blades, which are usually provided with finger and thumb ring handles. The distal end of the leaves can be designed differently. It is a reusable instrument intended for transient use.
Ear Loops; Nasal Poly-pus Snares; Tonsil Snares	A surgical ENT instrument, usually consisting of a canal or cannula and two fixed and sliding finger rings. A sling positions a circle of round, flexible wire around tissue in the ENT area and holds it in place. By moving the sliding finger ring you can control how much wire slides through the channel. The instrument is commonly used to exfoliate depleted foreign bodies. It is a reusable instrument intended for transient use.
Erisophakes	An ophthalmic product for rinsing the eye area and the eye socket without substantial change in energy, materials or parameters. It is a reusable, invasive instrument with respect to body orifices, which are not intended for connection to an active device and for transient use.
Scissors (tube shaft), endoscopy	A general, surgical instrument for cutting anatomical structures during a surgical procedure. It consists of ring handles, a tube shaft and a working end with fine cutting edges, e.g. straight, angled, or curved. It is made of stainless steel. It is a reusable instrument and intended for temporary use.
Punch (tube shaft)	A general, surgical instrument used to remove or remove tissue during a surgical procedure. It consists of ring handles, a tube shank and a working end with fine punches. It is made of stainless steel. It is a reusable instrument and intended for temporary use.
Pliers (tube shaft)	A general, surgical instrument for gripping, holding and manipulating anatomical structures during a surgical procedure. It consists of ring handles, a tubular shaft and a working end with fine jaws, e.g. straight, angled, or curved. It is made of stainless steel. It is a reusable instrument and intended for temporary use.

3.2 Indication

The use of tube shaft instruments is mainly in the following areas:

Grasping instruments	The instrument is used for treatments in gynecology. The procedure may only be carried out by suitably trained and qualified specialist personnel.
Erisophakes	The instrument is used for treatments in ophthalmology. It may only be used by trained and qualified professionals.
Prostatic Dilators	The instrument is used during examinations in urology. The procedure may only be carried out by suitably trained and qualified specialist personnel.
Haemorrhoidal Ligators; rectal snares	The instrument is used for patients with haemorrhoids. The treatment must be carried out by trained and qualified personnel.
Rectal Biopsy Specimen Forceps; Biopsy Punches	The instrument is used in various procedures for histopathological examinations. It is not intended for use in direct contact with the heart, central circulatory system, or central nervous system. The procedure may only be carried out by suitably trained and qualified specialist personnel.

<p>Bone Shears; Rib Shears</p>	<p>The instrument is used in various surgical procedures in orthopedics. The procedure may only be carried out by suitably trained and qualified specialist personnel.</p>
<p>Cotton Applicators; Needle holders with T/C jaws; Needle holders; Micro needle holders; Needle holders (tube shaft); Scissors (tube shaft), endoscopy; Punch (tube shaft); Pliers (tube shaft)</p>	<p>The instrument is used in various surgical procedures. The procedure may only be carried out by suitably trained and qualified specialist personnel.</p>
<p>Trocars; Breast Enlargement Sets; Ligature Guides; Ligature Conductors; Tendon passing instrument; Wire Strut Guide</p>	<p>The instrument is used in various surgical procedures. It is not intended for use in direct contact with the heart, central circulatory system, or central nervous system. The procedure may only be carried out by suitably trained and qualified specialist personnel.</p>
<p>Trachea Trocars; Scissors (Tube shaft); Goiter forceps; throat forceps; laryngoscopy instruments; laryngeal polyp forceps; nose and ear tweezers; nasal septum forceps; nose pliers; Ear forceps polyp forceps; Sharp nose pliers; septum forceps; Sphenoid Punches; Ethmoid Forceps; Antrum Punch; Nasal Scissors (Tube shaft); Septum Scissors; Nasal Scissors; Ear Loops; Nasal Polypus Snares; Tonsil Snares</p>	<p>The instrument is used during surgical ENT procedures. The procedure may only be carried out by suitably trained and qualified specialist personnel.</p>

3.3 Contraindications

No contraindications are known.

3.4 Patient target group

The products are suitable for all patient groups.

4 Warnings

	Medical products are delivered in a non-sterile condition and must be cleaned, disinfected and sterilized prior to their initial use.
	The use of faulty instruments is in principle forbidden and they have to go through the whole cleaning process before return.
	Please take into consideration that through higher power a bigger damage of the tissue can result: f.e. on forceps: the power at the end of the jaw is higher than at the tip of the jaw
	Remove all protective sleeves and films prior to first using or preparation for use.
	The safe combination of different products or of products with implants must be reviewed prior to clinical application by the user.
	Avoid improper throwing or dropping of instruments.
	Avoid mechanical overstressing of the instrument beyond the structural design, this can lead to breakage and deformation!
	A visual inspection of the instrument for damage and contamination must be carried out before each use!
	To prevent all contact corrosion, instruments with damaged surfaces must be separated immediately.
	If the products are used on patients with transmissible spongiform encephalopathy or HIV infection, we decline any responsibility for their reuse.
	After ophthalmical use, please pay attention to water quality during treatment (according to the specifications of AAMI TIR34 and the recommendations of the Rober-Koch-Institute on preparation of medical devices)!
	Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

5 Handling

The type of treatment must be determined in each individual case by the surgeon in cooperation with the internist and the anaesthetist.

For operational use in various surgical disciplines must be done by appropriately trained and qualified personnel.

6 Preparation

The person in charge of preparatory treatment is responsible for ensuring that the treatment is duly carried out using the relevant equipment, materials and personnel in the treatment facility and so achieves the desired result. This necessitates validation and routine monitoring of the process used. We urge you to take note of the national regulations dealing with instrument preparation.

The validated parameters refer to reusable surgical instruments. The validated parameters should also be observed for the other products described, unless a different procedure is explicitly described.

6.1 Information on instrument preparation

- Use cleaning and/or disinfection agents with a pH-value within 9-10.
- Please observe manufacturer instructions regarding dosage, exposure time and renewal of solutions.
- Use suitable brushes that do not remove material, i.e. **no** hard brushes (such as metal brushes and metal sponges) or coarse abrasive cleaners.
- Never leave instruments in cleaning or disinfection agents for longer than the specified time.
- Rinse and dry carefully through channels and pipes.
- Sensitive instruments must be cleaned in a storage or clamping fixture.
- Observe manufacturer instructions of cleaning – and sterilizing equipment.

6.1.1 Water quality

For cleaning, neutralization and rinsing, we recommend the use of demineralized water in accordance with the "Guideline DGKH, DGSV, AKI for the validation and routine monitoring of automated cleaning and thermal disinfection processes for medical devices and on the principles of device selection" (the guideline refers to DIN EN ISO 15883-1 point 6.4.2).

According to the guideline, the following values are recommended:

- Conductivity: $\leq 15 \mu\text{s}/\text{cm}$ (deviating from the table in DIN EN 285)
- pH value: 5 - 7
- Total hardness: $\leq 0.02 \text{ mmol CaO}/\text{l}$
- Salt content: $\leq 10 \text{ mg}/\text{l}$
- Phosphate (as P_2O_5): $\leq 0.5 \text{ mg}/\text{l}$
- Silicate (as SiO_2): $\leq 1 \text{ mg}/\text{l}$

- Chloride: ≤ 2 mg/l

6.2 Preparation at the place of use

Immediately after use, remove coarse dirt from the instruments with a suitable brush (see chapter 6.1) and rinse out the working cannulas. Do not use fixing agents or hot water ($> 40^{\circ}\text{C}$), as this results in residues becoming fixed and can affect the success of the subsequent cleaning operation

Dismantle and/or open instruments as far as possible. Within short time after use the instruments clean the instruments for reducing a drying of the residues.

This enables an easier cleaning. If instruments come into contact with corroding medicines or cleaning agents, wash these up with water immediately after use.

Longer drying times, e.g. for dry disposal are not validated and therefore not recommended.

V The drying time during validation was 1 hour.

6.3 Ultrasound bath (optional)

All instruments must be opened, dismantled and any cavities rinsed through.

Place instruments in the screen basket in such a way that overlaps and contact between instruments are avoided. Add cleaning agent to the water and adjust the temperature of the solution in line with the cleaning agent manufacturer's instructions.

The cleaning in the ultrasound bath should be at **35-40 kHz, 5 minutes** at least.

V To validate cleaning in an ultrasonic bath, the test items were ultrasonically treated in Neodisher Mediclean forte 0,5 % for 5 minutes.

Subsequently rinse instruments including all cavities before cleaning and disinfection.

The ultrasonic bath should not be used for medical devices where material could be removed in the ultrasonic bath, e.g. products made of soft plastic or coated instruments.

6.4 Manual cleaning

Since mechanical processes can be standardized, reproduced and therefore validated, mechanical cleaning/disinfection should be preferred to manual processes. Manual cleaning and disinfection process is not validated and therefore needs to be validated additionally by the end user.

6.5 Mechanical cleaning



Due to international standards (EN ISO 15883) and national guidelines, only validated mechanical cleaning and disinfection procedures should be used. For automated cleaning, we recommend a standard program for surgical instruments, e.g. instruments from Miele.

Observe the following by loading:

- Place the disassembled or opened instruments securely in the tray.
- Instruments with openings and recesses must be placed with the open side facing downwards so that they can be cleaned and no water from the cleaning process can collect in them.
- If available, use a coordinated rinsing device.

- Do not overload trays, avoid creating any overlaps


- Place the instruments with joints in an opened position into the cleaning, -and disinfection machine
- If a device for flushing, e.g. a Luer-Lock connection, is available, it should be used as

shown in the illustration..

- If needed use an adapter for the cleaning.

Preliminary rinsing is followed by chemical.

The chemical cleaning should take place at **40°C -60°C for at least 5 minutes.**

We recommend products with a **pH-value within 9-10**, e.g. Neodisher MediClean forte from Dr. Weigert. The cleaning agents used should be selected depending on the material and properties of the instruments and in accordance to national regulations: If there is a high chloride concentration in the water, pitting and tension crack corrosion can occur on the instruments.

The occurrence of this type of corrosion is minimized by using alkaline cleaning agents and demineralized water. By adding an acid-based neutralization medium, the rinsing off of alkaline cleaning agent residues is facilitated during the first intermediate rinsing process.

In order to prevent the formation of deposits, it is advisable to use neutral cleaners where the water quality is unfavourable.

After the second intermediate rinsing process, thermal disinfection takes place.

The thermal disinfection should be carried out with demineralized water that corresponds to microbiological drinking water quality (< 100 CFU/ml drinking water) at **80 to 95°C** and an exposure time in accordance with **EN ISO 15883**.

After the finishing of the programme take the good out of the machine because corrosion can arise if the instrument remains in the machine.

V Parameters used for the validation of preparation	
Pre-rinsing	1 minute with cold tap water
Cleaning	Temperature: 55 °C
	Soaking Time: 5 minutes (worst case)
	Neodischer Mediclean forte 0,4% (worst case)
Neutralization	Temperature: cold DI water
	Soaking Time: 2 minutes
	Neodisher Z 0,1%
Post-rinsing	2 minutes with cold DI water
Disinfection	Temperature: 90 °C (A ₀ 3000)
	Soaking Time: 5 minutes

6.6 Drying

Ensure adequate drying by the cleaning and disinfection device or using other suitable measures.

V Drying was omitted in the validation (worst case condition).

7 Maintenance, inspection

After cooling to room temperature, the instruments must be visually inspected for protein residues and other contamination. Slits, barriers, locks, tubes and other areas that are difficult to access must be thoroughly inspected. Instruments that are not residue-free must be repeatedly subjected to the entire reprocessing process.

To ensure that surgical instruments can be used for their intended purpose after reprocessing, it is necessary to carry out a functional test after cleaning, disinfection and drying of the visual inspection and care measures. Carry out the functional tests described in chapter 7.1.

Instruments that are stained, blunt, bent, no longer functional, damaged in any other way or the instrument labeling is no longer recognizable (e.g. labeling, UDI code no longer legible) must be discarded!

To help identify faulty instruments that need to be sorted out, we recommend the brochure "Instrument Reprocessing" from the Working Group "Instrumenten Aufbereitung". This includes Chapter 8 "Checks and Care" and Chapter 12 "Surface Changes: Deposits, Discoloration, Corrosion, Aging, Swelling and Stress Cracks".

7.1 Functional check

A newly purchased product must be subjected to a thorough visual and function check after its delivery and before each use.

Products must be checked for irregularities. Paying attention to cracks, fractures and the occurrence of corrosion.

If there are joints, the instruments should be oiled with a care product before the functional test. We recommend a medical white oil based on paraffin oil.

Check instruments with joints for ease of movement.

The labeling must also be checked for legibility.

Carry out a function check in accordance with the intended application of the instrument.

Essential tests for instruments with serration, jaw or cutting edges are among others:

- Correct opening and closing (smooth movement, complete)
- Intact serration (bent, broken off)
- The jaws must close properly
- No throwing the jaws over (cross bite)
- Smooth cutting edges
- No "scraping" or grinding noise when closing

Defective products must not be used and must have undergone the complete preparatory treatment process again before being returned.

8 Sterilization

Prior to sterilization, products must undergo cleaning and disinfection, be rinsed off without residue using demineralized water and subsequently dried. HEBUmedical recommends using a validated steam sterilization process (e.g. sterilizer in compliance with EN 285 and validated in accordance with DIN EN ISO 17665-1).

The validated parameters refer to reusable surgical instruments. The validated parameters should also be observed for the other products described, unless a different procedure is explicitly described.

On using the fractionated vacuum method, **sterilization** must be performed with at least **134°C (USA 132°C)** with a **minimum dwell period of 3 minutes**. Vacuum drying must then be carried out for at least 20 minutes.

V Parameters used for the validation of steam sterilization	
Prevacuum	3 times
Sterilization temperature	132 °C
Sterilization time	1,5 minutes (half cycle method)
Drying time	Drying time

The vapour must be free of ingredients, recommended limiting values of feed water and vapour condensate are determined through EN 285.

Other sterilization processes are compatible but not validated from HEBUmedical.

When loading, observe the recommended total weight. After the sterilization, check the sterile product packaging for damage, and inspect the sterilization indicators.

8.1 Packaging

Compliant packaging of products for sterilization in line with ISO 11607. Packaging used must be suitable for the instruments and protect them from microbiological contamination during storage. The seal must not be under tension. HEBUmedical recommends container or hospital common sterilization paper/film packagings for sterilized packaging.

V During validation the instruments were packaged in hospital common sterilization packagings (paper/film packagings) and steam sterilized.

9 Lifetime

The steam sterilization procedure was validated by laboratory tests. The products were sterile validated at a pre-vacuum of at least 5min duration and a temperature of 134°C for a lifetime of 50 cycles.

You can continue to use the instruments at your own responsibility over this cycle value if the tests described in chapter 7 have been successfully completed.

10 Storage

Store products in a dry, clean and dust-free environment at moderate temperatures from 5°C to 40°C.

Protect from the effects of the sun's rays and artificial light.

11 Warranty / Repair

Our products are manufactured from high-grade materials and carefully checked prior to dispatch. However, even if used properly in accordance with their intended purpose they are subject to a greater or lesser degree of wear depending on their intensity of use.

This wear is technically induced and unavoidable.

Should faults occur independently of wear, please contact our customer services. Defective products should no longer be used.

They must undergo the complete preparatory treatment process before being returned.

12 Disposal

Before disposal, the product must undergo the complete reprocessing and sterilization process in accordance with the procedure described in these instructions for use.

	<p>Risk of infection from contaminated products!</p> <ul style="list-style-type: none"> ■ When disposing of the product, its components and their packaging, the national regulations must be observed.
	<p>Risk of injury from sharp-edged and/or pointed products!</p> <ul style="list-style-type: none"> ■ When disposing of the product, ensure that the packaging prevents injury from the product.

13 Service and manufacturer address

Should you require the instructions for use in paper form, please use the contact details below. The instructions for use in paper form will be made available to you within seven calendar days of receipt of the request.

Alternatively, you can print out the electronic instructions for use yourself.

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HEBUmedical GmbH