Rohrschaftinstrumente / Tube shaft instruments

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Seite 5

Operating Manual

Page 15

DEUTSCH

ENGLISH

Inhalt

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Contents

1	Symbols	15
2	Introduction	15
3	Scope.....	15
3.1	Inteded Use	15
3.2	Contraindications.....	17
4	Warning.....	17
5	Handling.....	17
6	Preparation	17
6.1	Reutilization restrictions.....	18
6.2	Information on instrument preparation	18
6.3	Preparation at the place of use	18
6.4	Ultrasound bath (optional)	18
6.5	Manual cleaning	18
6.6	Mechanical cleaning	19
6.7	Drying.....	20
7	Sterilization	21
8	Maintenance, inspection.....	21
9	Packaging	21
10	Storage	21
11	Warranty / Repair	21
12	Manufacturer and servicing address.....	22

1 Symbols

Symbol	Definition
	CE marking
	Attention
	Manufacturer

2 Introduction

By purchasing this instrument, you are now the owner of a high-quality product whose use and correct handling are described in the following. In order to minimize possible risks to patients and users, please observe these instructions carefully. Use, disinfection, cleaning and sterilization may only be performed by suitably trained specialist personnel.

3 Scope

The use of tube shaft instruments takes place in almost every area of surgery, including:

- Abdominal surgery
- anesthesia
- ENT
- Hematology
- Urology

3.1 Inteded Use

Suction Tubes; rectal aspiration cannula: A unit of products for the periodic or continuous rinsing and / or aspiration of body cavities or wounds during a surgical procedure, without causing a significant change in energy, substances or parameters. It removes loose tissue, fluids, and other surgical field debris to keep it clear and improve vision and access. The unit usually includes a main powered pump, collection bottles, hoses and connections. It can be used with a rinsing catheter. It is not suitable for the supply of parenteral fluids, supply of medicines or other life-support measures. It is a reusable product.

Probes: A narrow, rod-like, flexible metal surgical instrument, usually with a blunt, beaded, or pointed tip, for examining fistulas, lumens, or wounds during an ENT procedure. The distal end, in front of the bead / tip, may be bent at right angles to the shaft. It is a reusable product.

Rectal director: A surgical instrument used for exploring the tissues and structures of the rectum during a gastroenterological / urological surgical procedure. This is a reusable instrument.

Ligature Conductors: A surgical instrument for guiding wire or ligatures through tissue. Typical guides may have flexible or rigid working parts with hooks, clips or claws for holding the material while passing through the tissue. The working end can have a button or an eye. Guides may also be instruments with ring handles, with straight or slightly curved working ends. The version with ring handles has a termination near the distal end. The working ends form short, grasping arms. It is a reusable instrument

Gall Stone Probes: A surgical instrument consisting of a slender rod of flexible metal with a blunt buttoned tip used to examine (bite) the bile duct, e.g. is used on the union of the ductus sympatheticus with the ductus cysticus. It is used during surgical procedures to examine or dilate bile duct narrowing. It is a reusable instrument.

Uterine Sounds: A surgical instrument for probing and measuring the internal length or depth of the uterus, cervix and vagina. It is a narrow, hollow or solid instrument made of steel or plastic in a cylindrical shape and is usually provided with length gradations on its working end. It is produced in different designs and elasticities. It is a reusable instrument.

Lachrymal Probes: A slender rod-shaped instrument made of flexible metal, with a blunt tip of tuber, used to examine eye tissue. It is a reusable product.

Cannula: A narrow, tapered metal tube for taking tissue samples for analysis. It is made in different diameters, the cannula tip is made in different versions (for example, straight or angled). It is a reusable product.

Erisophakes: An ophthalmic product for rinsing the eye area and the eye socket with a liquid to obtain a rinsing effect during an ophthalmic procedure and also to be able to aspirate fluid from the relevant area without substantial change in energy, materials or parameters. The product is used to remove debris, tissue or fluid from the surgical field, thereby helping to keep the surgical field clear and improve vision. Usually, the product is used in a cataract operation. The product works with a very low pressure to prevent accidental injury to the surgical field. It is a reusable product

Trocars: A surgical instrument with a sharp pyramidal or conical tip used to pierce body spaces. It is usually assembled and used with a compatible sleeve, thus allowing insertion of the instrument. After piercing, the trocar is withdrawn and provides a temporary working channel into the body cavity.

Trocars are usually rigid. It is a reusable product.

Trachea Trocars: A rigid surgical instrument with a sharp pyramidal or conical tip used to puncture body cavities during ear, nose and throat (ENT) surgery. It can be used in conjunction with a compatible sleeve into which it is inserted. Following puncture by the trocar located in the sleeve, the trocar is withdrawn, thereby releasing a temporary working channel through the sleeve into the body cavity. It is a reusable instrument.

Metal Catheters: A cylindrical rod-like surgical instrument for probing and / or stretching the urethra. It has a semi-rigid, wedge-shaped shaft that ends at the distal end in a narrow, beaded tip. It can be hollow to pass other narrow instruments. It is a reusable instrument

Rectal snares: A surgical instrument used to prevent hemorrhoids. A ligator usually consists of a long shaft or cannula with a clamp or handle at the distal end to hold the ligature. The proximal end consists of a handle with moving parts that mechanically actuate the distal clamp or handle. It is a reusable product.

Tendon passing instrument: A surgical instrument for guiding wire or ligatures through tissue. Typical guides may have flexible or rigid working parts with hooks, clips or claws for holding the material while passing through the tissue. The working end can have a button or an eye. Guides may also be instruments with ring handles, with straight or slightly curved working ends. The version with ring handles has a termination near the distal end. The working ends form short, grasping arms. It is a reusable instrument.

Rectal Biopsy Specimen Forceps: A surgical instrument for the general removal of biopsy samples from tumors and other tissue during surgical procedures for histopathological examinations. It is typically made of high quality stainless steel and is a long, inelastic instrument, the distal ends of which are provided with two shaped cutting shells with sharp edges, one exactly fitting into the other, or scissors-like shells with a cutting edge along its edge which will cut out the Allow biopsy sample. These are operated via scissor-like ring handles at the proximal end. It is introduced through an artificial or natural body opening in body cavities. The instrument is available in a variety of sizes. The instrument is reusable.

Cotton applicators: A product for use in combination with a suitable endoscope. It is used to bring a piece of absorbent gauze or fabric into the surgical site. It is a reusable product.

Ear Loops: A surgical ENT instrument, usually consisting of a canal or cannula and two fixed and sliding finger rings. A sling positions a circle of round, flexible wire around tissue in the ENT area and holds it in place. By moving the sliding finger ring you can control how much wire slides through the channel. The instrument is commonly used to exfoliate depleted foreign bodies. It is a reusable instrument

3.2 Contraindications

There are no known contraindications.

4 Warning

 <ul style="list-style-type: none"> Medical products are delivered in a non-sterile condition and must be cleaned, disinfected and sterilized prior to their initial use.
 <ul style="list-style-type: none"> The use of faulty instruments is in principle forbidden and they have to go through the whole cleaning process before return.
 <ul style="list-style-type: none"> Please take into consideration that through higher power a bigger damage of the tissue can result: f.e. on forceps: the power at the end of the jaw is higher than at the tip of the jaw
 <ul style="list-style-type: none"> Please observe the additional information enclosed with the products.
 <ul style="list-style-type: none"> Remove all protective sleeves and films prior to first using or preparation for use.
 <ul style="list-style-type: none"> The safe combination of different products or of products with implants must be reviewed prior to clinical application by the user.
 <ul style="list-style-type: none"> Avoid improper throwing or dropping of instruments
 <ul style="list-style-type: none"> To prevent all contact corrosion, instruments with damaged surfaces must be separated immediately.
 <ul style="list-style-type: none"> If used on patients with Creutzfeldt-Jakob disease or HIV infection, we are unable to accept any responsibility for reuse.

5 Handling

The type of treatment must be determined in each individual case by the surgeon in cooperation with the internist and the anaesthetist.

For operational use in various surgical disciplines must be done by appropriately trained and qualified personnel.

6 Preparation

The person in charge of preparatory treatment is responsible for ensuring that the treatment is duly carried out using the relevant equipment, materials and personnel in the treatment facility and so achieves the desired result. This necessitates validation and routine monitoring of the process used. We urge you to take note of the national regulations dealing with instrument preparation.

6.1 Reutilization restrictions

Frequent repeat preparatory treatment has minimal effects on the product. The end of the product life is normally determined by wear and damage due to use

6.2 Information on instrument preparation

- Use cleaning and/or disinfection agents with a pH-value within 9-10. Please observe manufacturer instructions regarding dosage, exposure time and renewal of solutions.
- Do not use hard brushes or coarse abrasive cleaners.
- Never leave instruments in cleaning or disinfection agents for longer than the specified time.
- Only used demineralized water for rinsing.
- Rinse and dry carefully through channels and pipes.
- Sensitive instruments must be cleaned in a storage or clamping fixture.
- Observe manufacturer instructions of cleaning – and sterilizing equipment.

6.3 Preparation at the place of use

Directly after using remove coarse dirt of the instruments and rinse out the working cannulas. Do not use fixing agents or hot water (> 40°C), as this results in residues becoming fixed and can affect the success of the subsequent cleaning operation

Dismantle and/or open instruments as far as possible. Within short time after use the instruments clean the instruments for reducing a drying of the residues. **Use cleaning and/or disinfection agents with a neutral pH-value <10.**

This enables an easier cleaning. If instruments come into contact with corroding medicines or cleaning agents, wash these up with water immediately after use.

6.4 Ultrasound bath (optional)

All instruments must be opened, dismantled and any cavities rinsed through.

Place instruments in the screen basket in such a way that overlaps and contact between instruments are avoided. Add cleaning

agent to the water and adjust the temperature of the solution in line with the cleaning agent manufacturer's instructions.

The cleaning in the ultrasound bath should be at 35-40 kHz, 5 minutes at least.

Subsequently rinse instruments including all cavities before cleaning and disinfection.

Medicine products which are possessing a bad ultrasound transmission, e.g. soft material are not usable for the ultrasound bath.

6.5 Manual cleaning

We recommend using lint-free soft cloths, plastic brushes and/or cleaning guns for manual cleaning. The cleaning agent should not have any fixing properties. The cleaning fluid must be changed daily.

Hints:

- If possible, apply the cleaning agent into cavities at pressure and carefully brush out.
- Be careful with the cleaning of barriers, slits, tight entries and joints etc.
- It has to be secured that instruments with joints will be cleaned in an opened and closed position.
- Use adequate brushes for the cleaning of cannulas, working channels and sack holes
- Flexible (complexe) instruments with unvisible surface and chuck have to be cleaned manually.
- We recommend for all push shafts, -and pipe shafts instruments and instruments whose surfaces are on top of each others during the cleaning (f.e.bone forceps and gouge forceps) a manually cleaning for an optimal cleaning result without residues

After that rinse the instruments with fully demineralized water extensively. Subsequently dip the instrument in disinfectant and use suitable cleaning instruments such as synthetic soft brushes. Observe the cleaning agent manufacturer's instructions. At the end of the cleaning process, repeat the rinsing action with fully demineralized water and dry the instruments completely.

6.6 Mechanical cleaning

On the basis of international standards (EN ISO 15883) and national directives, only validated machine cleaning and disinfection methods may be used. For the mechanical cleaning we recommend a standard programme for surgical instruments, f.e. instruments from Miele.

Observe the following by loading:



- Place the dismantled/opened instruments securely in the tray.
- Instruments with openings and gaps have to be faced down with the opened side so that they can be cleaned and no water of the cleaning process is collecting inside them. If available use balanced devices for rinsing

- Place the instruments with joints in an opened position into the cleaning,-and disinfection machine.
- If needed use an adapter for the cleaning

- Do not overload trays, avoid creating any overlaps.

Preliminary rinsing (cold, if applicable fully demineralized water without additives) is followed by chemical.

The chemical cleaning should take place at 40°C -60°C for at least 5 minutes.

We recommend products with a **pH-value within 9-10**, e.g. MediClean forte from Dr.Weigert. The cleaning agents used should be selected depending on the material and properties of the instruments and in accordance to national regulations: If there is a high chloride concentration in the water, pitting and tension crack corrosion can occur on the instruments. The occurrence of this type of corrosion is minimized by using alkaline cleaning agents and demineralized water. By adding an acid-based neutralization medium, the rinsing off of alkaline cleaning agent residues is facilitated during the first intermediate rinsing process (warm or cold water).

In order to prevent the formation of deposits, it is advisable to use neutral cleaners where the water quality is unfavourable. After the second intermediate rinsing process, thermal disinfection takes place.

The thermal disinfection should take place at temperatures of between 80 and 95°C, with an exposure time as outlined in EN ISO 15883.

After the finishing of the programme take the good out of the machine because corrosion can arise if the instrument remains in the machine.

6.7 Drying

Ensure adequate drying by the cleaning and disinfection device or using other suitable measures.

The drying temperature should not exceed 95°C to avoid material-related ageing processes.

7 Sterilization

Prior to sterilization, products must undergo cleaning and disinfection, be rinsed off without residue using demineralized water and subsequently dried. HEBUmedical recommends using a validated steam sterilization process (e.g. sterilizer in compliance with EN 285 and validated in accordance with DIN EN ISO 17665-1).

On using the fractionated vacuum method, **sterilization** must be performed with at least **134°C (USA 132° C)** with a **minimum dwell period of 4 minutes** and subsequent vacuum drying.

The vapour must be free of ingredients, recommended limiting values of feed water and vapour condensate are determined through EN 285.

Other sterilization processes are compatible but not validated from HEBUmedical.

When loading, observe the recommended total weight. After the sterilization, check the sterile product packaging for damage, and inspect the sterilization indicators

8 Maintenance, inspection

After cooling to room temperature, the instruments must be visually inspected for protein residues and other contamination, paying particular attention to cavities, blocks, inclusions, pipes, and other inaccessible areas. Instruments which are not free of residues must be returned for a complete retreatment process.

Carry out the functional check mentioned above. Instruments with stains, which are blunt, bent, no longer function or which are otherwise damaged must be segregated!

9 Packaging

Compliant packaging of products for sterilization in line with ISO 11607. Packaging used must be suitable for the instruments and protect them from microbiological contamination during storage. The seal must not be under tension. HEBUmedical recommends container for sterilized packaging.

10 Storage

Store products in a dry, clean and dust-free environment at moderate temperatures from 5°C to 40°C. Protect from the effects of the sun's rays and artificial light.

11 Warranty / Repair

Our products are manufactured from high-grade materials and carefully checked prior to dispatch. However, even if used properly in accordance with their intended purpose they are subject to a greater or lesser degree of wear depending on their intensity of use.

This wear is technically induced and unavoidable.

Should faults occur independently of wear, please contact our customer services. Defective products should no longer be used.

They must undergo the complete preparatory treatment process before being returned.

12 Manufacturer and servicing address

HEBUmedical GmbH
Badstraße 8
78532 Tuttlingen / Germany
Tel. +49 7461 94 71 - 0
Fax +49 7461 94 71 - 22
eMail: service@HEBUmedical.de
Web: www.HEBUmedical.de

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**Our vast range of products covers over 10 000 various instruments.
They can be found in Internet or request our catalog.**

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