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Contents

1	Symbols	13
2	Introduction	13
3	Scope	13
3.1	Inteded Use	13
3.2	Contraindications.....	14
4	Warning.....	14
5	Handling	15
6	Preparation.....	15
6.1	Reutilization restrictions	15
6.2	Information on instrument preparation	15
6.3	Preparation at the place of use	15
6.4	Ultrasound bath (optional)	16
6.5	Manual cleaning	16
6.6	Mechanical cleaning.....	16
6.7	Drying.....	18
7	Sterilization.....	18
8	Maintenance, inspection	18
9	Packaging	18
10	Storage.....	18
11	Warranty / Repair	19
12	Manufacturer and servicing address.....	19

1 Symbols

Symbol	Definition
	CE marking
	Attention
	Manufacturer

2 Introduction

By purchasing this instrument, you are now the owner of a high-quality product whose use and correct handling are described in the following. In order to minimize possible risks to patients and users, please observe these instructions carefully. Use, disinfection, cleaning and sterilization may only be performed by suitably trained specialist personnel.

3 Scope

The use of Standard Surgical Instruments is mainly in the following areas:

- abdominal surgery
- Spine surgery
 - Laminectomy
- thoracic surgery
 - Sternotomy
- Urology

3.1 Inteded Use

Anal Retractor: An instrument for temporarily holding back the walls of the anal canal. It is a reusable instrument.

Laminectomy spreaders, Lamina spreaders: A self-locking surgical instrument for temporarily separating the edges or withdrawing the wound edges of a surgical incision to permit the removal of a lamina of the vertebral arch. This may be necessary to remove a tumor, detachment, or tension caused by fracture or disc protrusion. The instrument resembles a large scissors with bulging branches with fenestrated, curved hooks that often end sharply to grab tissue. It is usually made of stainless steel and is available in different sizes and designs. It is a reusable instrument.

Cervical Traction Tongs: A surgical tool for temporarily gripping the skull to hold the head and neck in a specific, fixed position during a surgical procedure. It is a reusable instrument.

Retractor System: A frame-shaped device (e.g., a large ring) in which surgical sutures are hooked. The retractors may be dull or hooked, sometimes with several hooks. The retractors are commonly used to temporarily restrain the edges of a cut, or to temporarily hold off an organ or tissue during a surgical procedure. It is a reusable product.

Abdominal Retractor: A self-locking surgical instrument used to temporarily separate the wound edges during an abdominal cavity opening. It is commonly used during deep abdominal surgery, e.g. to make organs accessible or uncovered. The instrument may resemble a very large pair of scissors or be framed with blunt wide-angle blades which may be highly curved to form a saddle-shaped hook and / or wire-like blades that hold the abdomen open. It may have an additional third sheet stretched using a series of holes or notches. The instrument is made of high quality stainless steel. This is a reusable instrument.

Retractor: A locking device for temporary removal of tissue (e.g., keeping a wound open) or other anatomical structures or for exposing internal tissue or organs during surgery. The product resembles a pair of scissors with wide-angled leaves of various shapes (e.g., blunt, sharp, fenestrated, firm, deep). It is usually made of metal and reusable after sterilization.

Retractors; Amputating Retractor: A surgical instrument to temporarily separate and / or retract tissue to allow access to skeletal structures for examination or intervention. It is a reusable instrument.

Rib Contractors; Rib spreaders; rib spreaders made of Aluminium: A self-locking, frame-like surgical instrument for temporarily spreading two ribs to allow examination or intervention in the chest cavity. The instrument usually has two interchangeable leaves that hang over the ribs to be separated. Through a locking mechanism along a horizontal rail, the leaves can be detected and held in this position. The instrument is also known as a rib retractor and is manufactured in various designs and sizes. It is usually made of stainless steel, some variants for pediatrics are made of rust-resistant aluminum. It is a reusable instrument.

Bladder retractor: A frame-like mechanism to which the wound spreader is attached to temporarily hold the anatomical structures of the bladder. The wound spreaders are either blunt or hook-like (sometimes with several hooks). The wound spreader blades are used to temporarily restrain the edges of, for example: cuts, an organ, or tissue during a surgical procedure. It is a reusable instrument.

Endospecula: A self-locking surgical instrument designed to temporarily dislodge the vagina and temporarily retain the tissue to facilitate a surgical procedure or examination. The instrument may resemble a large pair of scissors or have an open circular-like shape, with tapered, blunted or wire-like leaves or grippers for grasping and dividing the tissue of the vaginal entrance. It is usually made of stainless steel and is available in different versions. It is a reusable instrument.

Lateral Vaginal Retractor: A product used during examination and treatment of the vagina to provide temporary access to the interior. Usually this is achieved by stretching the opening. The product is reusable and sterilizable.

Mouth Gags: An instrument that temporarily holds the mouth in an open position during an oral or dental procedure. It typically has blades and handles in various finishes, some models have additional features to support stability. It is a reusable instrument.

Nasal Specula: A surgical instrument, usually composed of two articulated branches, available in a variety of sizes, shapes and contours. It has a blunt distal end and is used in insertion and opening to stretch or stretch tissue in the nostril. It is used for examination or during procedures. It is usually made of stainless steel with a self-opening mechanism or an adjustable mechanism.

3.2 Contraindications

There are no known contraindications.

4 Warning

Medical products are delivered in a non-sterile condition and must be cleaned, disinfected and sterilized prior to their initial use.

The use of faulty instruments is in principle forbidden and they have to go through the whole cleaning process before return.

Please take into consideration that through higher power a bigger damage of the tissue can result: f.e. on forceps: the power at the end of the jaw is higher than at the tip of the jaw

!	Please observe the additional information enclosed with the products.
!	Remove all protective sleeves and films prior to first using or preparation for use.
!	The safe combination of different products or of products with implants must be reviewed prior to clinical application by the user.
!	Avoid improper throwing or dropping of instruments
!	To prevent all contact corrosion, instruments with damaged surfaces must be separated immediately.
!	If used on patients with Creutzfeldt-Jakob disease or HIV infection, we are unable to accept any responsibility for reuse.

5 Handling

The type of treatment must be determined in each individual case by the surgeon in cooperation with the internist and the anaesthetist.

For operational use in various surgical disciplines must be done by appropriately trained and qualified personnel.

6 Preparation

The person in charge of preparatory treatment is responsible for ensuring that the treatment is duly carried out using the relevant equipment, materials and personnel in the treatment facility and so achieves the desired result. This necessitates validation and routine monitoring of the process used. We urge you to take note of the national regulations dealing with instrument preparation.

6.1 Reutilization restrictions

Frequent repeat preparatory treatment has minimal effects on the product. The end of the product life is normally determined by wear and damage due to use

6.2 Information on instrument preparation

- Use cleaning and/or disinfection agents with a pH-value within 9-10. Please observe manufacturer instructions regarding dosage, exposure time and renewal of solutions.
- Do not use hard brushes or coarse abrasive cleaners.
- Never leave instruments in cleaning or disinfection agents for longer than the specified time.
- Only used demineralized water for rinsing.
- Rinse and dry carefully through channels and pipes.
- Sensitive instruments must be cleaned in a storage or clamping fixture.
- Observe manufacturer instructions of cleaning – and sterilizing equipment.

6.3 Preparation at the place of use

Directly after using remove coarse dirt of the instruments and rinse out the working cannulas. Do not use fixing agents or hot water (> 40°C), as this results in residues becoming fixed and can affect the success of the subsequent cleaning operation

Dismantle and/or open instruments as far as possible. Within short time after use the instruments clean the instruments for reducing a drying of the residues. **Use cleaning and/or disinfection agents with a neutral pH-value <10.**

This enables an easier cleaning. If instruments come into contact with corroding medicines or cleaning agents, wash these up with water immediately after use.

6.4 Ultrasound bath (optional)

All instruments must be opened, dismantled and any cavities rinsed through.

Place instruments in the screen basket in such a way that overlaps and contact between instruments are avoided. Add cleaning

agent to the water and adjust the temperature of the solution in line with the cleaning agent manufacturer's instructions.

The cleaning in the ultrasound bath should be at 35-40 kHz, 5 minutes at least.

Subsequently rinse instruments including all cavities before cleaning and disinfection.

Medicine products which are possessing a bad ultrasound transmission, e.g. soft material are not usable for the ultrasound bath.

6.5 Manual cleaning

We recommend using lint-free soft cloths, plastic brushes and/or cleaning guns for manual cleaning. The cleaning agent should not have any fixing properties. The cleaning fluid must be changed daily.

Hints:

- If possible, apply the cleaning agent into cavities at pressure and carefully brush out.
- Be careful with the cleaning of barriers, slits, tight entries and joints etc.
- It has to be secured that instruments with joints will be cleaned in an opened and closed position.
- Use adequate brushes for the cleaning of cannulas, working channels and sack holes
- Flexible (complexe) instruments with unvisible surface and chuck have to be cleaned manually.
- We recommend for all push shafts, -and pipe shafts instruments and instruments whose surfaces are on top of each others during the cleaning (f.e.bone forceps and gouge forceps) a manually cleaning for an optimal cleaning result without residues

After that rinse the instruments with fully demineralized water extensively. Subsequently dip the instrument in disinfectant and use suitable cleaning instruments such as synthetic soft brushes.


Observe the cleaning agent manufacturer's instructions. At the end of the cleaning process, repeat the rinsing action with fully demineralized water and dry the instruments completely.

6.6 Mechanical cleaning



On the basis of international standards (EN ISO 15883) and national directives, only validated machine cleaning and disinfection methods may be used. For the mechanical cleaning we recommend a standard programme for surgical instruments, f.e. instruments from Miele.

Observe the following by loading:

- Place the dismantled/opened instruments securely in the tray.
- Instruments with openings and gaps have to be faced down with the opened side so that they can be cleaned and no water of the cleaning process is collecting inside them. If available use balanced devices for rinsing

- Place the instruments with joints in an opened position into the cleaning,-and disinfection machine.
- If needed use an adapter for the cleaning
- Do not overload trays, avoid creating any overlaps.

Preliminary rinsing (cold, if applicable fully demineralized water without additives) is followed by chemical.

The chemical cleaning should take place at 40°C -60°C for at least 5 minutes.

We recommend products with a **pH-value within 9-10**, e.g. MediClean forte from Dr.Weigert. The cleaning agents used should be selected depending on the material and properties of the instruments and in accordance to national regulations: If there is a high chloride concentration in the water, pitting and tension crack corrosion can occur on the instruments. The occurrence of this type of corrosion is minimized by using alkaline cleaning agents and demineralized water. By adding an acid-based neutralization medium, the rinsing off of alkaline cleaning agent residues is facilitated during the first intermediate rinsing process (warm or cold water).

In order to prevent the formation of deposits, it is advisable to use neutral cleaners where the water quality is unfavourable. After the second intermediate rinsing process, thermal disinfection takes place.

The thermal disinfection should take place at temperatures of between 80 and 95°C, with an exposure time as outlined in EN ISO 15883.

After the finishing of the programme take the good out of the machine because corrosion can arise if the instrument remains in the machine.

6.7 Drying

Ensure adequate drying by the cleaning and disinfection device or using other suitable measures.

The drying temperature should not exceed 95°C to avoid material-related ageing processes.

7 Sterilization

Prior to sterilization, products must undergo cleaning and disinfection, be rinsed off without residue using demineralized water and subsequently dried. HEBUmedical recommends using a validated steam sterilization process (e.g. sterilizer in compliance with EN 285 and validated in accordance with DIN EN ISO 17665-1).

On using the fractionated vacuum method, **sterilization** must be performed with at least **134°C (USA 132° C)** with a **minimum dwell period of 4 minutes** and subsequent vacuum drying.

The vapour must be free of ingredients, recommended limiting values of feed water and vapour condensate are determined through EN 285.

Other sterilization processes are compatible but not validated from HEBUmedical.

When loading, observe the recommended total weight. After the sterilization, check the sterile product packaging for damage, and inspect the sterilization indicators

8 Maintenance, inspection

After cooling to room temperature, the instruments must be visually inspected for protein residues and other contamination, paying particular attention to cavities, blocks, inclusions, pipes, and other inaccessible areas. Instruments which are not free of residues must be returned for a complete retreatment process.

Carry out the functional check mentioned above. Instruments with stains, which are blunt, bent, no longer function or which are otherwise damaged must be segregated!

9 Packaging

Compliant packaging of products for sterilization in line with ISO 11607. Packaging used must be suitable for the instruments and protect them from microbiological contamination during storage. The seal must not be under tension. HEBUmedical recommends container for sterilized packaging.

10 Storage

Store products in a dry, clean and dust-free environment at moderate temperatures from 5°C to 40°C. Protect from the effects of the sun's rays and artificial light.

11 Warranty / Repair

Our products are manufactured from high-grade materials and carefully checked prior to dispatch. However, even if used properly in accordance with their intended purpose they are subject to a greater or lesser degree of wear depending on their intensity of use.

This wear is technically induced and unavoidable.

Should faults occur independently of wear, please contact our customer services. Defective products should no longer be used.

They must undergo the complete preparatory treatment process before being returned.

12 Manufacturer and servicing address

HEBUmedical GmbH
Badstraße 8
78532 Tuttlingen / Germany
Tel. +49 7461 94 71 - 0
Fax +49 7461 94 71 - 22
eMail: service@HEBUmedical.de
Web: www.HEBUmedical.de

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**Our vast range of products covers over 10 000 various instruments.
They can be found in Internet or request our catalog.**

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